

# **Alaris™ Infusion Central v1.2**

## **Quick Start Guide**

14 May 2019

**ascom**



*This guide is provided with the product Alaris Infusion Central to give the user essential information on the following topics:*

- *Alaris Infusion Central start up and user access*
- *Frequently used procedures*
- *Recovery procedures*
- *User manual availability*
- *Product information availability*
- *Manufacturer and Distributors contacts*

# Contents


<b>Contents .....</b>	<b>3</b>
<b>“Alaris Infusion Central” installation.....</b>	<b>4</b>
<b>Frequently used procedures .....</b>	<b>4</b>
How to run “Alaris Infusion Central” .....	4
Log In procedure .....	5
Module Selection .....	6
Log out.....	6
Patient Management .....	7
Infusion module .....	10
Devices management.....	16
Fluid Balance module .....	18
System procedures.....	19
<b>System unavailability .....</b>	<b>20</b>
<b>“Alaris Infusion Central” user manual .....</b>	<b>20</b>
On line update of the user manual .....	21
Printed user manual .....	21
<b>Contacts .....</b>	<b>23</b>
BD contacts .....	23
Manufacturer .....	24


## Alaris Infusion Central installation

The Alaris Infusion Central system can only be installed by Ascom Ums/Distributors technicians or technicians explicitly authorized by Ascom UMS or Distributors. The end-user will find the system already installed and properly configured on his/her workstation.

## Frequently used procedures

### How to run Alaris Infusion Central

To run Alaris Infusion Central Double click the  icon on the workstation desktop. The following window is displayed while the system is loading.



After a few seconds the Alaris Infusion Central work area is displayed.

AIC	Choose patient...	ICU CENTRAL	08:48 18 MAR 2019
1 PATIENT01, TEST01 01	2 PATIENT02, TEST02 02 0h 12m to next end of infusion	3 PATIENT03, TEST03 03 0h 12m to next end of infusion	4 PATIENT04, TEST04 04 0h 12m to next end of infusion
Alaris GP 2.25 mL/h	Alaris GP 2.25 mL/h	Alaris GP 2.25 mL/h	Alaris GP 2.25 mL/h
Alaris GW 5 mL/h	Alaris GW 5 mL/h	Alaris GW 5 mL/h	Alaris GW 5 mL/h
Alaris GH_G 30 mL/h	Alaris GH_G 30 mL/h	Alaris GH_G 30 mL/h	Alaris GH_G 30 mL/h
Alaris GH_G 9 mL/h	Alaris GH_G 9 mL/h	Alaris GH_G 9 mL/h	Alaris GH_G 9 mL/h
Alaris CC 0.9 mL/h	Alaris CC 0.9 mL/h	Alaris CC 0.9 mL/h	Alaris CC 0.9 mL/h
Alfentanil 11.2 mL/h	Alfentanil 11.2 mL/h	Alfentanil 11.3 mL/h	Alfentanil 11.2 mL/h
5 PATIENT05, TEST05 05 0h 12m to next end of infusion	6 PATIENT06, TEST06 06 0h 12m to next end of infusion	7 PATIENT07, TEST07 07 0h 12m to next end of infusion	8 PATIENT08, TEST08 08 0h 07m to next end of infusion
Alaris GP 2.25 mL/h	Alaris GP 2.25 mL/h	Alaris GP 2.25 mL/h	Alaris GP 3.6 mL/h
Alaris GW 5 mL/h	Alaris GW 5 mL/h	Alaris GW 5 mL/h	Alaris GW 5 mL/h
Alaris GH_G 30 mL/h	Alaris GH_G 30 mL/h	Alaris GH_G 30 mL/h	Alaris GH_G 30 mL/h
Alaris GH_G 9 mL/h	Alaris GH_G 9 mL/h	Alaris GH_G 9 mL/h	Alaris GH_G 9 mL/h
Alaris CC 0.9 mL/h	Alaris CC 0.9 mL/h	Alaris CC 0.9 mL/h	Alaris CC 0.9 mL/h
Alfentanil 11.3 mL/h	Alfentanil 11.2 mL/h	Alfentanil 11.2 mL/h	Alfentanil 9.27 mL/h
PROFILE	DOSE	RATE	PRESSURE
VOLUME	TIME	WEIGHT	ROTATE

## Log In procedure

To use Alaris Infusion Central it is necessary to log in to the system. I.e. the user must insert his/her credentials (USERNAME and PASSWORD). To log in, at the beginning of every work session,

1. Click the **User** button (A). The following screen is displayed.

The screenshot shows the login interface of Alaris Infusion Central. At the top, a blue command bar contains a 'User' button (A), a 'Choose patient...' button (B), and a 'ICU CENTRAL' button (C). Below the command bar, there are input fields for 'USERNAME' and 'PASSWORD'. A numeric keypad is displayed below the password field. At the bottom, there is a 'RECENT' list with a table of recent users. The 'OK' button (D) is highlighted in the bottom right corner.

RECENT				
1	ADMIN	2		3
4		5		6
7		8		9
10		11		12

2. Insert your **USERNAME** (B).
3. Insert your **PASSWORD** (C).
4. Click **OK** (D).

The user is logged in. An acronym indicating the user currently logged is displayed on the **User** button on the command bar.



## “User lock” functionality

If the system remains idle for a certain time (configurable), the user is automatically logged out. To disable this functionality it is necessary, before clicking **OK**, to:

- Click **LOCK** on the “Log in” screen command bar (**A**).



If the user is locked, a padlock icon is displayed on the User icon.



## Module Selection

To select a module click the corresponding icon on the lateral bar.

The selected module's icon is highlighted (**A**).



## Log out

At the end of every work session the user must log out. To log out:

- Click the **User** button on the command bar (**A**).



The user is logged out. The user acronym disappears from the **User** button -



## Patient Management

### ◆ Patient selection

1. Click the **Patient** button (A) on Control Bar.



2. Select **PATIENTS ADMITTED** (B) and then **IN BED** (C)



3. Click the **Bed Button** corresponding to the patient to be selected (D).




The patient is selected. Their name is displayed on the **Patient** button on “Control Bar”. Click **None** on the command bar (E) to deselect the selected patient.


### ◆ New patient creation and admission

If the system is connected to the hospital patients' archive, the patient can be admitted and selected using the procedure described in the next paragraph. If the patient cannot be found in the hospital patients' archive, it is possible to create it in the Alaris Infusion Central local database using the procedure described in this paragraph.

1. Click **NEW/ADMIT** on the command bar (A).



The following window opens.



**NEW/ADMIT PATIENT**

Family Name:  Given Name:  Initials:

Patient Code:  Birth Date:  Sex:

Notes:

Admission Date - time: 25/01/2019 13:54:35 Discharge Date - time:

Admission Code:  Height [cm]:  Weight [kg]:

Location: ICU Bed: 6

OK CANCEL

Callout B points to the Patient Code field. Callout C points to the Location and Bed fields.


- Specify the patient data (B).
- Specify the bed and location (C).
- Click **OK**.

#### ◆ Patient admission from the hospital patient archive

- Click the **Patient** button (A) on “Control Bar”.



- Select **PATIENTS SEARCH** (B) and then **REMOTE SEARCH** (C)



**PATIENTS SEARCH**

**REMOTE SEARCH**

First name:  Last name: JOHN Birth date: / / Sex: M Patient code:

SEARCH CLEAR

First name	Last name	Sex	Birth date	Patient code
DOE	JOHN	M	23/06/1998	333333
SMITH	JOHN	M	12/12/1977	111111

LOCAL SEARCH REMOTE SEARCH ANONYMOUS

Callout B points to the PATIENTS SEARCH button. Callout C points to the REMOTE SEARCH button.

- Specify patient data in the search fields (D) and click **SEARCH** (E).



- Double-click the row corresponding to the relevant result (F).

First name: DOE Last name: J Birth date: / / Sex: F Patient code:

SEARCH CLEAR

First name	Last name	Sex	Birth date	Code
DOE	JANE	F	12/12/1912	564738

- Verify patient data. Specify the destination department and bed (G) for the patient and click **OK**.

**NEW/ADMIT PATIENT**

Family Name: DOE Given Name: JOHN Initials:

Patient Code: 45376 Birth Date: 12/08/1977 Sex:

Notes:

Admission Date - time: 18/03/2019 09:04:05 Discharge Date - time:

Admission Code: Height [cm]: Weight [kg]:

Location: ICU Bed: 8

OK CANCEL

### ◆ Patient discharge

Select the patient as described in the previous paragraphs.

- Click **EDIT (A)** on the command bar. A menu is displayed.

PATIENT EXPLORER + NEW/ADMIT EDIT NONE CLOSE

- On the menu, click the **DISCHARGE (B)** option.

EDIT

DELETE

DISCHARGE

ADMIT

MOVE

A confirmation is requested. Click **Yes** to confirm. A window containing the patient data opens.

EDIT PATIENT

Family Name

Stanley

Given Name

Initials

Patient Code

57663

Birth Date

09/09/2018

Sex

M

Notes

Admission Date - time

11/09/2018 11:58:18

Discharge Date - time

28/01/2019 08:53:22

Admission Code

Height [cm]

186

Weight [kg]

81.000

OK

CANCEL

3. Specify the discharge data and click **OK (C)**.

## “Infusion” module

### ◆ Ward station

The “Ward station” screen displays all the pumps connected to each patient in the domain.

<b>1 PATIENT01, TEST01 01</b> 	<b>2 PATIENT02, TEST02 02</b> 0h 54mto next end of infusion <div> <div>Alaris GP</div> <div>3.6 mL/h</div> </div> <div> <div>Alaris GW</div> <div>7 mL/h</div> </div> <div> <div>Alaris GH_G</div> <div>30 mL/h</div> </div> <div> <div>Alaris GH_G</div> <div>9 mL/h</div> </div> <div> <div>Alaris CC</div> <div>1.2 mL/h</div> </div> <div> <div>Alfentanil</div> <div>6.64 mL/h</div> </div>	<b>3 PATIENT03, TEST03 03</b> 0h 54mto next end of infusion <div> <div>Alaris GP</div> <div>3.6 mL/h</div> </div> <div> <div>Alaris GW</div> <div>7 mL/h</div> </div> <div> <div>Alaris GH_G</div> <div>30 mL/h</div> </div> <div> <div>Alaris GH_G</div> <div>9 mL/h</div> </div> <div> <div>Alaris CC</div> <div>1.2 mL/h</div> </div> <div> <div>Alfentanil</div> <div>6.64 mL/h</div> </div>	<b>4 PATIENT04, TEST04 04</b> 0h 54mto next end of infusion <div> <div>Alaris GP</div> <div>3.6 mL/h</div> </div> <div> <div>Alaris GW</div> <div>7 mL/h</div> </div> <div> <div>Alaris GH_G</div> <div>30 mL/h</div> </div> <div> <div>Alaris GH_G</div> <div>9 mL/h</div> </div> <div> <div>Alaris CC</div> <div>1.2 mL/h</div> </div> <div> <div>Alfentanil</div> <div>6.64 mL/h</div> </div>
<b>5 PATIENT05, TEST05 05</b> 0h 54mto next end of infusion <div> <div>Alaris GP</div> <div>3.6 mL/h</div> </div> <div> <div>Alaris GW</div> <div>7 mL/h</div> </div> <div> <div>Alaris GH_G</div> <div>30 mL/h</div> </div> <div> <div>Alaris GH_G</div> <div>9 mL/h</div> </div> <div> <div>Alaris CC</div> <div>1.2 mL/h</div> </div> <div> <div>Alfentanil</div> <div>6.64 mL/h</div> </div>	<b>6 PATIENT06, TEST06 06</b> 0h 54mto next end of infusion <div> <div>Alaris GP</div> <div>3.6 mL/h</div> </div> <div> <div>Alaris GW</div> <div>7 mL/h</div> </div> <div> <div>Alaris GH_G</div> <div>30 mL/h</div> </div> <div> <div>Alaris GH_G</div> <div>9 mL/h</div> </div> <div> <div>Alaris CC</div> <div>1.2 mL/h</div> </div> <div> <div>Alfentanil</div> <div>6.64 mL/h</div> </div>	<b>7 PATIENT07, TEST07 07</b> 0h 54mto next end of infusion <div> <div>Alaris GP</div> <div>3.6 mL/h</div> </div> <div> <div>Alaris GW</div> <div>7 mL/h</div> </div> <div> <div>Alaris GH_G</div> <div>30 mL/h</div> </div> <div> <div>Alaris GH_G</div> <div>9 mL/h</div> </div> <div> <div>Alaris CC</div> <div>1.2 mL/h</div> </div> <div> <div>Alfentanil</div> <div>6.64 mL/h</div> </div>	<b>8 PATIENT08, TEST08 08</b> 0h 49mto next end of infusion <div> <div>Alaris GP</div> <div>3.6 mL/h</div> </div> <div> <div>Alaris GW</div> <div>7 mL/h</div> </div> <div> <div>Alaris GH_G</div> <div>30 mL/h</div> </div> <div> <div>Alaris GH_G</div> <div>9 mL/h</div> </div> <div> <div>Alaris CC</div> <div>1.2 mL/h</div> </div> <div> <div>Alfentanil</div> <div>5.96 mL/h</div> </div>

PROFILE

DOSE

RATE

PRESSURE

VOLUME

TIME

WEIGHT

ROTATE

The screen is divided in rectangular areas. Every area, called “bed-card”, refers to a bed and contains a schematic representation of all the pumps connected to a bed.

On top of each bed-card the bed number and patient name are displayed (A). Below the patient name the time remaining until the next end of infusion is specified (B). The rows indicated in C represent the connected pumps. Each row represents a pump. The rows can be displayed in five colors:

Bed Number	Patient Name	Time to Next End of Infusion	Pump Name	Rate (mL/h)
1	John, Doe JD	20m to next end of infusion	Amiodaron	20
			Noradrenaline	23
			Alaris CC	45
			Frusemide	22
			Midazolam	10
			Alaris GW	30
			DOPAMINE	0
			DOPAMINE	10


#### ◆ “Infusion” module color code

On the “Infusion” module the following color code applies:

State	Example	Icon
Stand-by	Pump is in stand-by.	
Infusing	Pump is infusing.	
Low priority alarm	Pump is sending a low-priority alarm.	
Medium priority alarm	Pump is sending a medium-priority alarm.	
High priority alarm	Pump is sending a high-priority alarm.	




## ◆ “Infusion” module icons

Different icons can be displayed on the “Bed cards” (A).



1. Click the icons to open a tooltip providing additional information.

The icons meaning and number is set in the configuration. Contact your system administrator for more information. Alaris Infusion Central has three pre-configured icons.

-  - *Patient weight mismatch/profile mismatch.* This icon is displayed either when different weights for the same patient are set on the pumps in which the dose rate is influenced by the patient weight, or when different profiles are specified on different pumps for the same patient.
-  - *Expired vascular access device.* This icon is displayed when one of the vascular access devices associated to the patient exceeds the maximum number of days set on the Invasive Device Management module.
-  - *Guided Relay.* This icon is displayed either when one or more of the administered drugs are eligible for a Guided Relay process or a Guided Relay is running.

Different icons are displayed on the left of each pump name (B) to identify different kinds of pumps (volumetric, syringe, enteral, target-controlled infusion etc...).

## ◆ Patient Station

Click on one of the “Bed cards” to open the “Patient station” screen. The “Patient station” screen offers a detailed view of all the data coming from the pumps connected to a patient. The corresponding patient is automatically selected. On the left of the screen is a list of syringes and infusion pumps connected to the patient (A); in the middle a diagram displays drug infusion velocity changes in time and possible administered boluses (B).



On the left, each box represents a pump. These boxes are named “Pump buttons”. The pump button displays the drug name when the pump provides this kind of information. When it doesn’t, the pump name is displayed.

A notification area is displayed on the right of every Alaris Infusion Central screen, reporting the possible alarms sent by the connected pumps (C).


The notification area can be, by configuration:

- Always visible
- Automatically displayed when a new notification arrives
- Only visible after user click on the **Notification** button on “Control Bar”.

The different notifications are displayed in chronological order, (most recent on top) and by criticality.

◆ Pump detail

On the “Patient station” screen, click one of the pump buttons to display a screen containing detailed information on the pump.



On the left, a list of all the events occurred on the selected pump is displayed (A). On the right, three charts are displayed, representing some of the trends of the current infusion parameters: rate, dose, volume, pressure (B).


## ◆ Guided Relay Process

The “Guided Relay” feature is designed to support clinicians in transitioning continuous infusion of critical drugs which come to end and need to be changed with as minor an impact as possible in blood’s drug concentration.

This task is accomplished by pre-setting a gradual change of syringe, with one syringe gradually diminishing the infusion dose rate and another syringe (infusing the same drug at the same time) gradually increasing the infusion dose rate.

The dose rate changes (defined as “steps”) must be performed on the pump by the nursing staff. Alaris Infusion Central provides a series of reminders that can be of aid in deciding the time and amount of the dose rate changes in the two syringes. If the Alaris Infusion Central notifications are ignored, the Guided Relay process is automatically aborted.

When the Guided Relay process is set up, the infusion charts change according to the specified Guided Relay settings (see figure below).



Each step of the process is displayed on the chart, on the right of the yellow now-bar. The required changes to the dose rate are indicated for each step (**A**). See the Alaris Infusion Central user manual for the Guided Relay setup procedure.



## Devices management

### ◆ Add a new device

On the “Invasive Device Management” module, to add a new device

1. Under DEVICE, click **Add** on the command bar (A).



The following screen is displayed.

A screenshot of the 'ADD NEW DEVICE' form in the Alaris Infusion Central software. The form contains several input fields: 'Device:' (a dropdown menu with a plus icon), 'Site:' (a dropdown menu with a plus icon), 'Max days:' (a text input field), 'Activity:' (a dropdown menu with 'Insertion' selected), 'Date Time:' (a date and time picker showing '21/01/2019' and '13:50'), 'Operator:' (a dropdown menu with 'ADMIN' selected and a plus icon), and 'Notes:' (a large text area). At the bottom of the form, there are three buttons: 'Remove' (with a minus icon), 'Ok' (with a checkmark icon), and 'Cancel' (with an 'X' icon). A red callout bubble with the letter 'B' points to the 'Device:' dropdown menu. Another red callout bubble with the letter 'C' points to the 'Ok' button.

2. Specify the data of the new device (B).
3. Click **Ok** (C).



## ◆ Add a new nursing activity

1. Click, on the devices table of the “Invasive Device Management” module, the row corresponding to the device to which the nursing activity refers (A).

JOHN

Code: 1111, Birthdate: 07/09/1989, Age: 29 years, Sex: male  
Admission date: 11/09/2018, Days from admission: 132

LIST OF VASCULAR ACCESS DEVICES

Device (site)	Days (max)	Date Time (elapsed)	Operator	Activity (notes)
Device 1 (Site 1)	0 (3)	2019-01-21 01:50 (0 h)	ADMIN	Insertion (Notes on device)
Device 2 (Site 2)	0 (1)	2019-01-21 01:52 (0 h)	ADMIN	Present at admission
Device 3 (Site 3)	1 (3)	2019-01-20 12:55 (1 h)	NURSE 1	Insertion (Note)

COMMAND BAR: + Add | Edit | X Remove | + Add | Edit | Filter: All | InSitu | Rem... | Print | Expand

2. Under NURSING, click **Add** on the command bar (B).

The following screen is displayed.

ADD NEW NURSING

Device: Device 2

Site: Site 2

Max days: 1

Activity: Activity x

Date Time: 21/01/2019 13:58

Operator: ADMIN

Notes:


COMMAND BAR: X Remove | Ok | Cancel

3. specify the data of the new nursing activity (C).
4. Click **Ok** (D).

## “Fluid Balance” module

### ◆ Entering data in the fluid balance

1. Select the “Fluid Balance” module.
2. Click **New** on the command bar (A).



The following screen is displayed.

3. Select the Input/Output parameter (B).
4. Enter the parameter value (C).
5. Repeat step 3 and 4 for all the values to be entered.
6. Click **Save** to confirm (D).

## System procedures

### ◆ Quit Alaris Infusion Central

To exit from Alaris Infusion Central

1. click the **Menu** button on “Control bar” (A). The Alaris Infusion Central main menu opens.



2. Click **QUIT** on the menu. Another menu opens.
3. Click **QUIT AIC** to exit. Click **SHUT DOWN AND RESTART** to exit from Alaris Infusion Central and restart the workstation.


### ◆ Check sound and brightness

To check the functionality and level of sound and brightness during the night and day use the following procedure.

1. click the **Menu** button on “Control bar” (A). The Alaris Infusion Central main menu opens.



2. Click **SYSTEM CONFIGURATION**. Another menu opens.
3. Click **NIGHT/DAY TESTER**. A specific test window opens.
4. Select, on the drop-down menu indicated in the figure (B), the kind of notification to be tested.
5. Click **Check Night Mode** and **Check Day Mode** (C) to check sound and brightness in both modes.



In case brightness and/or sound are inappropriate please contact the system administrators to correct them.

## System unavailability

If, during start up, there are problems connecting with the server, the system provides a specific “System recovery” window.

The connection problem sometimes is automatically solved after a short time. If not, please contact Ascom UMS/Distributors technical assistance. See page 23 for Ascom UMS/Distributors contacts.

Cases exist, rare but possible, in which it is physically impossible to use the Alaris Infusion Central system.

The healthcare organization using Alaris Infusion Central should define an emergency procedure to face this kind of events and to:

- 1) allow the hospital departments to be operative again;
- 2) restore the system availability as soon as possible.


Ascom UMS can provide support and assistance for the emergency procedure, but the healthcare organization is responsible for the procedure definition and execution.

The contacts list is on page 23.

## Alaris Infusion Central user manual

The Alaris Infusion Central system user manual is provided with the Product. It can be displayed in digital format within the Alaris Infusion Central application. To access the user manual:

- Click the **Help** button on the command bar (A).



A screen opens, making it possible to access the Alaris Infusion Central user manuals. The on line documentation is in PDF format.

## On line update of the user manual

The Alaris Infusion Central user manual can also be downloaded from the internet. Every Alaris Infusion Central customer can access a protected area on which all the existing versions of the manuals are present. The user can this way check if an updated version of their manual exists.

To download the Alaris Infusion Central manuals, type the following address on your web browser address bar:

<https://confluence.ascom-ws.com/display/DIG/Alaris+Infusion+Central>

A welcome page is displayed. Enter username and password to access the User Manuals list.

**Username** – User name provided by Ascom UMS/Distributor

**Password** – Password provided by Ascom UMS /Distributor

## Printed user manual

The customer can request a printed version of the user manuals of the purchased Products. Ascom UMS/Distributors will provide the requested manuals as soon as possible. The request must be forwarded to the technical assistance. See “Contacts” section.

To display the Alaris Infusion Central about box

- 
 AIC
 Choose patient...
  ADM
  ICU CENTRAL
 08:56 18 MAR 2019
 

-  **BD Alaris™ Infusion Central 1.2** 
- Product:** DIGITAT<sup>®</sup>  
**Configuration:** Alaris Infusion Central  
**Version:** 6.0.0.4  
**Date:** 2019-05-15
-  *Instructions for use are supplied in electronic form. Requirements for viewing the instructions: PDF reader.  
Contact your BD representative for details.*
- TEMPORARY LICENCE - SERVER**  
Issued to: UMS on 05/22/19 Port: 50001  
Expiration date: 12/31/2019  
Message Dispatcher on localhost:50500  
SN: AICENT\_20190522\_5HWFE5A3HE4FP2499B655UX94KA665
-  **Manufactured by ASCOM UMS srl**  
**Via Amikare Ponchielli, 29**  
**50018 Scandicci (FI), Italy**
- Digitat<sup>®</sup> is a registered trademark of ASCOM UMS srl.  
Copyright (c) by ASCOM UMS srl. All rights Reserved.  
Portions Copyright © GrapeCity, Inc. 1987-2017. All Rights Reserved.  
Portions Copyright © Innu System 1994-2017. All Rights Reserved.  
Portions Copyright © Teleoh 2017. All Rights Reserved.*

## Contacts

### BD contacts

Country	Telephone	email
Australia	Freephone: 1 800 656 100	<a href="mailto:AUS_customerservice@bd.com">AUS_customerservice@bd.com</a>
België	+32(0)53 720 556	<a href="mailto:info.benelux@bd.com">info.benelux@bd.com</a>
Canada	+1 800 387 8309	<a href="mailto:CanadaCapital@carefusion.com">CanadaCapital@carefusion.com</a> .
Danmark	+45 43 43 45 66	<a href="mailto:bddenmark@bd.com">bddenmark@bd.com</a>
Deutschland	+49 622 1305 558	<a href="mailto:GMB-DE-CustService-Alaris@bd.com">GMB-DE-CustService-Alaris@bd.com</a>
España	+34 902 555 660	<a href="mailto:info.spain@carefusion.com">info.spain@carefusion.com</a>
France	+33 (0) 1 30 02 81 41	<a href="mailto:mms_infusion@bd.com">mms_infusion@bd.com</a>
Italia	+39 02 48 24 01	<a href="mailto:customer.service-italy@bd.com">customer.service-italy@bd.com</a>
Magyarország	(36) 1 488 0233	<a href="mailto:info.cfn.export@bd.com">info.cfn.export@bd.com</a>
Nederland	+31(0)20 582 94 20	<a href="mailto:info.benelux@bd.com">info.benelux@bd.com</a>
New Zealand	Freephone: 0800 572 468	<a href="mailto:NZ_customerservice@bd.com">NZ_customerservice@bd.com</a>
Norge	+47 64 00 99 00	<a href="mailto:bdnorge@bd.com">bdnorge@bd.com</a>
Polska	+48 22 377 11 00	<a href="mailto:Info_Poland@bd.com">Info_Poland@bd.com</a>
Portugal	+351 219 152 593	<a href="mailto:dl-pt-geral@carefusion.com">dl-pt-geral@carefusion.com</a>
South Africa	Freephone: 0860 597 572 +27 11 603 2620	<a href="mailto:bdsa_customer_centre@bd.com">bdsa_customer_centre@bd.com</a>
Suomi	+358-9-8870 780	<a href="mailto:bdsuomi@bd.com">bdsuomi@bd.com</a>
Sverige	+46 8-7755160	<a href="mailto:bdsweden@bd.com">bdsweden@bd.com</a>
Switzerland	+41 61 485 22 66	<a href="mailto:Customer_Care_Switzerland@bd.com">Customer_Care_Switzerland@bd.com</a>
United Kingdom	Freephone: 0800 917 8776	<a href="mailto:BDUK_CustomerService@bd.com">BDUK_CustomerService@bd.com</a>
United States of America	Freephone: 800 482 4822	<a href="mailto:CustCareInfusion@carefusion.com">CustCareInfusion@carefusion.com</a>
中国	400 878 8885	<a href="mailto:serviceclientbdf@bd.com">serviceclientbdf@bd.com</a>

**Manufacturer**

Ascom UMS srl unipersonale  
Via A. Ponchielli 29  
50018 Scandicci (FI)  
Italy  
+39 055 0512161  
[it.info@ascom.com](mailto:it.info@ascom.com)