



# **Smart Monitor Web User Manual**

Version 7.0

2023-03-30

# Contents

|  |           |
|--|-----------|
| <b>1. Smart Monitor Web.....</b>                 | <b>3</b>  |
| 1.1 Introduction.....                            | 3         |
| 1.2 Intended use .....                           | 3         |
| <b>2. Wearables.....</b>                         | <b>4</b>  |
| <b>3. Login .....</b>                            | <b>4</b>  |
| <b>4. Global Dashboard.....</b>                  | <b>5</b>  |
| 4.1 Creating a search filter.....                | 8         |
| <b>5. Patient Dashboard .....</b>                | <b>9</b>  |
| <b>6. Associating a kit to a patient.....</b>    | <b>11</b> |
| <b>7. Detaching a device from a patient.....</b> | <b>15</b> |
| <b>8. Association check.....</b>                 | <b>16</b> |
| <b>9. Wearable Devices .....</b>                 | <b>17</b> |
| <b>10. User-Location Configuration.....</b>      | <b>19</b> |
| <b>11. Patient admission mode .....</b>          | <b>20</b> |
| 11.1 Patient admission .....                     | 21        |
| 11.2 Patient Discharge .....                     | 22        |

## 1. Smart Monitor Web



For general and detailed information about the Product environment and the intended use of the product, see the specific documents. The knowledge and understanding of these documents is mandatory for an appropriate and safe use of the Smart Monitor Web, described in this document.

---



Smart Monitor Web will not be used for patient monitoring that could require immediate actions or direct clinical assistance. Smart Monitor Web shall be used for the clinical conditions of patients affected by chronic diseases that require a daily or periodic monitoring of some vital parameters.

---



The caregivers who are monitoring patient data should always verify and validate the measures (for ex. by requesting to repeat them, etc.) in order to minimize the risk of other family members and relatives using the same wearable device.

---

### 1.1 Introduction

Smart Monitor Web is a web application used to view data collected from medical devices, both wearables or devices directly connected to the beds of one or more patients, admitted to hospital facilities or remotely monitored for home therapy. The module provides a variously configurable dashboard that offers an overview of tens or hundreds of patients by displaying, for each patient, the latest available vital signs, and historical data.

Additional administrative functionalities are also available like the association between Gateway App and patients, the admission and discharge of patients.



*In some configurations Smart Monitor Web is directly connected to the medical devices (i.e. Patient Monitor). In those cases, the Gateway App is not present and the Patient-Kit association is not required.*

---

### 1.2 Intended use

Smart Monitor Web is part of the Digistat Suite. Digistat Suite is split in two different products according to the functionalities implemented in the different modules (see the product manuals for a more detailed description). Smart Monitor Web is part of Digistat Care and inherits the same intended use.

Read the product manuals (USR ENG Digistat Care) for a detailed description of the intended use and disclaimer notes.

The Digistat Suite can transfer the collected data to third party systems. Please, verify on the third-party system how the information provided by the Digistat Suite will be used.

## 2. Wearables

In addition to data from medical devices associated with hospital beds, the Smart Monitor Web can also display historical data coming from a set of “kits” connected to the patients. Every kit collects the data from a single patient.

A “kit” has the following components:

- Supported Android Smartphone (ex. Ascom Myco 3).
- Ascom Gateway App (see user manual USR ENG Gateway).
- One or more wearable devices (see above mentioned manual for the complete and updated list of supported wearable devices).

Once a kit has been assembled (a smartphone with the Gateway app installed and wearables) and configured, it is possible to associate it to a patient via Smart Monitor Web. After a kit is associated, it is possible to see the status of the collection of parameters in the main dashboard or in the patient detail dashboard.

Smart Monitor Web makes it also possible to disassociate the kit from the patient when data collection is no longer required.



---



The devices supported by Smart Monitor Web are: VivaLnk (VV330 and VV200), Biovotion Everion, iHealth (PO3, BP5, BP5s and BP550BT), Oxitone 1000M, Temp Sitter WT1, Gemini BP, Vitalograph Model 4000, VivaLnk Checkme SpO2 sensor, Cosinuss c-med alpha.

The units of measure supported by Smart Monitor Web are: “bpm” for Heart Rate; “r/min” for Respiratory Rate; “°C” for Temperature; “%” for SpO2; “mmHg” for Pressure, “Kg” for weight.

---

## 3. Login

To access the Smart Monitor Web module either select:

- this icon  from the Digistat ControlBar Desktop or this icon  from ControlBar Web;
- alternatively, from the Configurator Web > Web Modules page, copy on a browser or click on the **configured URL for Smart Monitor Web** for the browser navigation.

Perform the login entering the correct **username** and **password** and then, click **Login** to authenticate the user.

If an error is displayed, try again checking for typos. If the problem persists, contact system administrators. If authentication is successful, a Dashboard will be displayed to navigate the application.

## 4. Global Dashboard

The global dashboard displays the latest collected vital parameters for a set of patients connected to medical devices at the hospital bed or to wearable devices for home monitoring and admitted to configured locations. A logged user can display only the patients admitted to his “visible” locations (see chapter 10).

Drag a column header (Location) and drop it here to group by that column

| Location | Slot | Patient   | HR (bpm)  | RR (bpm)     | Skin Temp. (°C) | SpO2 (%) | NBPd (mm Hg) | NBPs (mm Hg) | Priority | INFO |
|----------|------|---|-----------|--------------|-----------------|----------|--------------|--------------|----------|------|
| ICU      | 2    | Cartonio Tonio<br>TC002 Age: 33 Years<br>Sex: ♂   | 63<br>now | 21.49<br>now | 89<br>now       | N/A      | 77<br>now    | 134<br>now   | N/A      | N/A  |
| ICU      | 3    | Cotogno Milo<br>MC003 Age: 23 Years<br>Sex: ♂     | 61<br>now | 21.02<br>now | 97<br>now       | N/A      | 83<br>now    | 125<br>now   | N/A      | N/A  |
| ICU      | 4    | Lina Fata<br>FL004 Age: 34 Years<br>Sex: ♀        | 70<br>now | 19.97<br>now | 24<br>now       | N/A      | 87<br>now    | 123<br>now   | N/A      | N/A  |
| ICU      | 5    | Lucente Amina<br>AL005 Age: 22 Years<br>Sex: ♀    | 63<br>now | 21.71<br>now | 87<br>now       | N/A      | 87<br>now    | 115<br>now   | N/A      | N/A  |
| ICU      | 6    | Rosarospa Strega<br>SR006 Age: 21 Years<br>Sex: ♀ | 71<br>now | 23.07<br>now | 95<br>now       | N/A      | 87<br>now    | 131<br>now   | N/A      | N/A  |
| ICU      | 8    | Manno Orco<br>OM008 Age: 36 Years<br>Sex: ♂       | N/A       | 23.76<br>now | N/A             | N/A      | N/A          | N/A          | N/A      | N/A  |
| ICU      | 1    | Fosco Lupo<br>LF001 Age: 40 Years<br>Sex: ♂       | 70<br>now | 20.32<br>now | 87<br>now       | N/A      | 85<br>now    | 119<br>now   | N/A      | N/A  |

7 rows

SMART MONITOR WEB

Admit Patient Discharge Patient

Fig 1

The main table displays the following information:

- **Location:** hospital/clinic/ward to which the patient is associated or admitted to (patients can be either physically located in a hospital or at home).
- **Slot:** it represents a bed, or a virtual position (slot) used to virtually locate the patient.
- **Patient:** patient personal data, as name and surname, ID, age, and sex.
- **List of vital sign parameters:** a configurable list of columns displaying some vital signs parameters, with different icons and colors. Every vital parameter is displayed with the following information: value and collection time (ex. 1 min ago, 1 hour ago). The columns can be customized through the **GridConfiguration** system option. See the *CFG ENG Digistat Suite MDR* manual for details on the columns' configuration for the main Dashboard.

The table can be sorted clicking on any header. Location columns support grouping by dragging the column header to the area on top of the table.




Fig 2

It is also possible to filter the table using the buttons located on the upper-left corner (Fig 1). So:

- Click on a provided button to filter the table according to a “customized” rule (Fig 2) associated to the button.




Fig 3

Filter buttons can be added, edited, or deleted, using the **Actions** menu in the upper right corner of the table (Fig 1 B), then clicking on the **Manage filter** button (Fig 3).

It is possible to search for a specific patient using the search filter (Fig 1 C) located in the upper part of the table:

- insert patient last name or first name, in full or even partially,
- then click on the **Search** button.

Fig 4

On the bottom part of the screen, a command bar is provided with different possible buttons enabled according to the configuration of **NoDeviceAssignment** system option (Fig 4). If **NoDeviceAssignment** is set to true, the device assignment won't be available and the buttons for patient admission and discharge will be provided as the example shows.

SMART MONITOR WEB

Drag a column header (Location) and drop it here to group by that column

| Location | Slot | Patient  | HR (bpm)  | RR (bpm)     | Skin Temp. (°C) | SpO2 (%) | NBPd (mm Hg) | NBPs (mm Hg) | Priority | INFO |
|----------|------|--|-----------|--------------|-----------------|----------|--------------|--------------|----------|------|
| ICU      | 2    | <b>Cartonio Tonio</b><br>TC002 Age: 33 Years<br>Sex: ♂   | 64<br>now | 19.6<br>now  | 36<br>now       | N/A      | 74<br>now    | 121<br>now   | N/A      | N/A  |
| ICU      | 3    | <b>Cotogno Milo</b><br>MC003 Age: 23 Years<br>Sex: ♂     | 61<br>now | 20.18<br>now | 40<br>now       | N/A      | 77<br>now    | 137<br>now   | N/A      | N/A  |
| ICU      | 4    | <b>Lina Fata</b><br>FL004 Age: 34 Years<br>Sex: ♀        | 65<br>now | 23.08<br>now | 37<br>now       | N/A      | 75<br>now    | 131<br>now   | N/A      | N/A  |
| ICU      | 5    | <b>Lucente Amina</b><br>AL005 Age: 22 Years<br>Sex: ♀    | 70<br>now | 20.35<br>now | 38<br>now       | N/A      | 85<br>now    | 128<br>now   | N/A      | N/A  |
| ICU      | 6    | <b>Rosarospa Strega</b><br>SR006 Age: 22 Years<br>Sex: ♀ | 68<br>now | 19.85<br>now | 39<br>now       | N/A      | 74<br>now    | 121<br>now   | N/A      | N/A  |
| ICU      | 7    | <b>Ronfo Gnomo</b><br>RO007 Age: 22 Years<br>Sex: ♂      | 70<br>now | N/A<br>now   | 38<br>now       | N/A      | 74<br>now    | 137<br>now   | N/A      | N/A  |

15 rows

Assign Device Detach Device Devices

Fig 5

Otherwise, if the system option value is not flagged and it is set to false, three buttons are provided instead of the “**Admit Patient**” and “**Discharge Patient**” ones (Fig 1 D):

- **Assign Device** (Fig 5 A): to start the device assignment procedure.
- **Detach Device** (Fig 5 B): to start the device detachment procedure.
- **Devices** (Fig 5 C): to read the list of devices available for the assignment and those already attached to a bed/patient.




Fig 6

For further information on the configuration of the command bar, read the document *DSO ENG System Options*.

Clicking on the three dots button, available in both scenarios, a menu is provided with two options:

- **Verify device association** (Fig 6 A): to check the associated patient.
- **Configure User-Location** (Fig 6 B): to create, edit, delete, or just view available users and their associated locations.

In the following paragraphs, the actions and functions mentioned up to now will be reviewed and explained in detail.

## 4.1 Creating a search filter

To create or edit the pre-defined search filters:

- click on the **Actions** (Fig 1 B) button, and
- select the **Manage Filters** (Fig 3) option.

The **Manage filters** window appears:




Fig 7

The existing filters are listed on a grid the left part of the window (Fig 7 A) and each filter has an index number. The order of the filters can be edited, selecting one at a time and then moving it up and down through the list using the **up and down buttons** arrows provided (Fig 7 B). The filters can also be deleted, selecting them and the red x button (Fig 7 C) provided for each of them. They can be edited, modifying the fields available on the right side of the window (**Error! Reference source not found. E**).

It is also possible to create a new filter and to do that:

- Click the **+** button (**Error! Reference source not found. D**) provided on top of the grid.
- Customize the filter filling in the available fields.
- Press **Save** to save the filter and close the filters window.

A filter has the following properties:

- **Name**: the name of the filter. This text is displayed on the button in the main dashboard. Consequently, it is better to enter short names.
- **Index**: it is the sorting index, not editable.



Fig 8

- **Expression:** it is possible to create a logical expression to specify the filter rules. It is possible to add as many logical constructs as needed (Fig 8 A). For example: the Test Filter displayed, once saved, and selected, will show only the patients with HR greater or equal to 100 bpm **and** temperature greater or equal to 38 C° **and** SpO2 less or equal to 93%. A **textual summary** is provided under the expression fields.

## 5. Patient Dashboard

Double clicking on a patient card, the user can access the **Patient Detail Dashboard**. Two possible scenarios are given:

Fig 9

- If **Online Web** has been installed, a page has been correctly configured and the ID of that page has been indicated in the **OnlinePageID** system option (Fig 9) of the WearableMonitor module, the patient detail dashboard coincides with that Online Web page and the user is redirected to the Online page. See *CFG ENG Online Validation* manual and the document *DSO ENG System Options* for details.

- If the **OnlinePageID** system option is set to 0 (the default value), the Smart Monitor Web patient detail page is given.



Fig 10

The patient detail dashboard page consists of several part that contains different information:

- A central section (Fig 10 A), that contains the **charts** returning and displaying the **trends** and **history** of the **vital parameters** retrieved over time for the selected patient. The charts display data in a 6-hour period by default. The x-axis represents the time while the y-axis represents the values. Click any chart to display the value at a certain time in a tooltip.
- On the top part (Fig 10 B), some buttons are given to change and customize the **time interval**. Use the buttons "6 hours", "12 hours" and "24 hours" to change the displayed time interval. It is also possible to manually specify a time value (number of hours) in the field provided and then press the **Search** button to apply the customization.
- On the right part of the page (Fig 10 C), all the **events** and **alarms** are listed, both those communicated by the patient (using the Gateway App) and those auto-generated (for example: "device disconnected").
- On the lower part, the **command bar** shows the **Dashboard button** (Fig 10 D), to go back to the main Dashboard view.
- Above the event table **four buttons** (Fig 10 E), are provided:



Fig 11


The grey buttons are disabled while the black ones are enabled. These buttons activate the following functionalities (left to right):

- **Tracking** (Fig 11 A): it displays on a map where the kit consisting of wearables devices, smartphone and Gateway app is located. To enable this functionality the Gateway app must be configured to read GPS coordinates.
- **Info** (Fig 11 B): it displays the assigned medical device data like the name and code/unique identifier.
- **Telephone** (Fig 11 C): it starts a telephone call with a configured phone number.
- **Patient Privacy Document** (Fig 11 D): it displays the **Patient Privacy Document**, generated during the association procedure (see chapter 6).

## 6. Associating a kit to a patient

In case the Gateway-Wearable kit combination is used, to start the patient-device association process, the smartphone must have been connected to the system at least once, in order to have it correctly registered. To start the association workflow:

- Click **Assign Device** button on Smart Monitor Web command bar (Fig 5 A). The following window is displayed:



The 'Assign Device' window is a web-based form with a black title bar and a close button. It contains four numbered steps:

- 1 SELECT PATIENT**: Includes a text input field for 'Insert Patient Code:' and a 'SEARCH' button.
- 2 SELECT DEVICE**: Includes a text input field for 'Insert Device Code:' and a 'SEARCH' button.
- 3 SELECT LOCATION**: Includes two dropdown menus, 'Select Location:' (with 'ICU' selected) and 'Select Bed:'.
- 4 PATIENT SIGNATURE**: Includes the instruction 'Please print the 'Delivery and Privacy Agreement' document and ask the patient to sign it!' and a 'PRINT DOCUMENT' button.

At the bottom right of the form area are two buttons: 'Assign device' and 'Close'.

Fig 12

The association workflow can be completed in four steps:

The 'Assign Device' window displays a four-step process. Step 1, 'SELECT PATIENT', is highlighted with a red 'A' icon. It includes a text input field containing 'AL005' and a 'SEARCH' button. To the right, a 'SELECTED PATIENT' box (marked with a red 'B' icon) shows patient details: 'LUCENTE, AMINA AL', 'Born on 12/04/2000 in N/A', and 'Code: AL005'. Step 2, 'SELECT DEVICE', has an empty 'Insert Device Code' field and a 'SEARCH' button. Step 3, 'SELECT LOCATION', features dropdown menus for 'Select Location' (set to 'ICU') and 'Select Bed' (set to '5'). Step 4, 'PATIENT SIGNATURE', includes a 'PRINT DOCUMENT' button and instructions to print a 'Delivery and Privacy Agreement' document. At the bottom right are 'Assign device' and 'Close' buttons.

Fig 13

### Step 1

Enter the **patient code** used in this specific installation and press the **SEARCH** button (Fig 13 **A**). If the patient already exists, then the related record is displayed in the **Selected Patient box** (Fig 13 **B**).

The 'Assign Device' window shows Step 1, 'SELECT PATIENT', with a red 'A' icon. The 'Insert Patient Code' field contains '123' and the 'SEARCH' button is visible. To the right, the 'PATIENT SEARCH' box (marked with a red 'A' icon) displays the message 'No patient found with code 123.' in red text.

Fig 14

If the patient code is not retrieved, a red written message (Fig 14 **A**) appears warning the user that no patient has been found.

**Assign Device**

**1 SELECT PATIENT**  
Insert Patient Code:  
AL005  
SEARCH

**2 SELECT DEVICE**  
Insert Device Code:  
SK00004841  
SEARCH

**3 SELECT LOCATION**  
Select Location: ICU  
Select Bed: 5

**4 PATIENT SIGNATURE**  
Please print the 'Delivery and Privacy Agreement' document and ask the patient to sign it!  
PRINT DOCUMENT

**SELECTED PATIENT**  
LUCENTE, AMINA AL ♂  
Born on 12/04/2000 in N/A  
Code: AL005

**SELECTED DEVICE**  
Gateway  
Code: SK00004841  
Status: AVAILABLE

Assign device Close

Fig 15

## Step 2

Enter the ID of the gateway app (the kit ID) or the ID of the medical device to be associated (Fig 15 A). Click on **SEARCH** button and if the selected device is retrieved, its data are displayed in a new box. The status is set to **Available** (Fig 15 B).



*The kit ID can be read on a label attached to the smartphone (if present) or in the Gateway app, selecting the nurse modality (see the relevant documentation for more instructions).*

If the ID is not found, then it is necessary to check for typos. If no typos are detected, in case the ID of a kit is searched, it may be the case that the kit is new to the network. In this case, connect the kit to the network and try again.

**2 SELECT DEVICE**  
Insert Device Code:  
SK00000266  
SEARCH

**FOUND DEVICE**  
Gateway  
Code: SK00000266  
Status: ASSIGNED TO PATIENT: Rossi, Mario - Code: 666999666  
DETACH DEVICE

Fig 16

If the device or kit are already associated to a patient, a red written message appears in the Status field (Fig 16), warning the user that they have been already assigned to a different patient, reporting their name, surname, and code. In this case, either quit the current workflow or press **DETACH DEVICE** to start the detachment workflow for the current device/kit. The detachment of the device from the original patient enables the subsequent assignment to the current patient.

## Step 3

If a patient is not already admitted to a bed, it is possible to select the location where the user wants the patient to be admitted and one free slot. It is possible to choose only from those locations that are enabled for current user.

If the selected patient was already admitted to a bed, this section is disabled, and user can proceed to step 4.

Fig 17

#### Step 4


This step is not mandatory but configurable through the

**MandatorySignatureOnDeviceAssignment** system option (Fig 17). If the value of the system option is set to true, during the device assignment procedure, the signature will be mandatory in order to proceed with the association, otherwise it won't be mandatory and the user will be able to associate a device without displaying the document.

Fig 18

So, if the **MandatorySignatureOnDeviceAssignment** system option is set to true, as in the example provided, to complete the association workflow (and to activate the ASSOCIATE DEVICE button) it is necessary to press the **PRINT DOCUMENT** button (Fig 15 C). A privacy report is displayed. If the specific healthcare organization procedures require it, it is possible to print the document for patient signature. After the document has been displayed, the **Assign Device** button enables. Click on it to confirm the association. See the document *DSO ENG System Options* for details on the configuration of the system option.


## 7. Detaching a device from a patient



The screenshot shows a window titled "Detach Device" with a close button (X) in the top right corner. Inside the window, there is a section labeled "SELECT DEVICE". Below this label, it says "Insert Device Code:". There is a text input field containing the code "IX3M3342". Below the input field is a button labeled "SEARCH". At the bottom of the window, there are two buttons: "Detach device and close" and "Close".

Fig 19

The procedure to detach a Gateway-Wearable kits or other medical devices from patients can be performed by clicking on **Detach Device** (Fig 5 B) button on Smart Monitor Web command bar. The **Detach Device** window is displayed:



The screenshot shows the same "Detach Device" window. The text input field now contains the code "IP4M13302". To the right of the input field, there is a red error message that says "No assigned device found with code IP4M13302." Two red callout boxes with yellow letters "A" and "B" are present. Callout A points to the text input field, and callout B points to the error message. The "SEARCH" button is still visible below the input field. The "Detach device and close" and "Close" buttons are at the bottom.

Fig 20

- Enter the **device ID** in the field provided (Fig 20 A) and press the **SEARCH** button. If the device is not found an error message (Fig 20 B) is displayed warning the user that no assigned device with the code entered has been retrieved, otherwise, the following window is displayed:




Fig 21

On the right section, a box that contains the **device information** (name and code, Fig 21 A) and the patient - to which it is assigned - **personal data** is displayed (Fig 21 B). Furthermore, a **checkbox** (Fig 21 C) is provided to discharge the patient at the same time as the finalization of the device detachment procedure.

- Press **Detach device and close** button to complete the disassociation workflow.

The device is then available for another patient.

If that device is the only one associated to the previous patient, then the patient data disappears from the main dashboard.

## 8. Association check

The **Check device-patient association functionality** makes it possible to verify if a device is correctly associated to a patient. This workflow performs a double check in order to reduce risks of a wrong association. It can be performed, for example, in case the association workflow is performed before the actual "live" association to a patient.




Fig 22

To check the association:

- Click on three dots button (Fig 22 A) on Smart Monitor Web command bar and then on **Verify device association** (Fig 22 B). The following window opens:






Fig 23

Type the device code in the field provided (Fig 23 **Error! Reference source not found. A**) and click on **Check associated patient** (Fig 23 **Error! Reference source not found. B**).

If an Android device (e.g. a smartphone) is used, it is possible to tap the **BARCODE scanner** button (Fig 23 **Error! Reference source not found. C**) to read the device barcode. To read the barcode it is necessary to install, on the smartphone, the "Barcode Scanner" third party app (by ZXing Team).




Fig 24

If the device is found and it is already assigned to a patient, then the associated patient data is displayed (Fig 24 **A**), otherwise an error message appears.

## 9. Wearable Devices

It is possible to display the list of all medical devices that have been connected to the system. To do that, just:



Fig 25

- Click on **Devices button** on Smart Monitor Web command bar (Fig 25 **A**). The following screen is displayed:

| Export to Excel  |                         |   |                         |                     |  |
|------------------|-------------------------|---|-------------------------|---------------------|--|
| Device Code      | Device Label            | Patient Family Name                               | Patient Name            | Last Connection     |  |
| IX665013         | test label test         | ***   | ***                     |                     |  |
| IX3M5013         | Test 2                  | Family20  | Given20                 | 07/07/2022 13:38:02 |  |
| IX6M5013         |                         | Simpson   | Richard                 | 07/07/2022 13:38:02 |  |
| IX9M5013         |                         | TEST  | PATIENT                 | 07/07/2022 13:38:02 |  |
| IX12M5013        |                         |   |                         |                     |  |
| IX15M5013        |                         | John  | Doe                     | 07/07/2022 13:38:02 |  |
| 12345            |                         | PatientNr6  | 5                       | 07/07/2022 13:38:02 |  |
| xxx              |                         |   |                         |                     |  |
| SK00003032       |                         | COGNOME628  | NAME628                 |                     |  |
| IX18M5013        |                         | ***   | ***                     | 07/07/2022 13:38:02 |  |
| IX21M5013        |                         | Dfgfg   | High2                   | 07/07/2022 13:38:02 |  |
| 720a7fe77d254f19 |                         | Patient   | 1                       |                     |  |
| IX24M5013        |                         | Palmer  | Laura                   | 07/07/2022 13:38:02 |  |
| IP1M12401        | DEV_LABEL2_NICOLA       | IamAVeryLongFamilyName<br>MadeByTwoDifferentParts | AlsotheGivenName IsLong | 01/10/2021 12:15:00 |  |
| IP1M12802        | IP1M12802               | John  | Doe                     | 07/07/2022 13:38:02 |  |
| IP2M12701        | Test                    | ***   | ***                     | 07/07/2022 13:38:02 |  |
| IP2M12402        | XXX678                  | Jordan  | Gregory                 | 07/07/2022 13:38:02 |  |
| IP2M13303        |                         | ***   | ***                     |                     |  |
| IP2M12404        | Test label associated 3 | White   | Emily                   |                     |  |
| IP2M12605        |                         | ***   | ***                     | 07/07/2022 13:38:02 |  |

Fig 26

The table displays the list of all the medical devices and Gateway-Wearables kits that, at least once, have been connected to the system.

The “**Device Code**” column contains the devices IDs. If a device is already associated to a patient, the “**Patient Family Name**” and “**Patient Name**” columns contain respectively the name and surname of the patient that is currently associated. **Device Label** and **Last Connection** columns are also provided, to respectively shows information on the Label used to customize the device with a simpler and more recognizable name, and, in case the device is currently in use, the date and time information about the last connection retrieved.

| Device Details   |           |                     |      |
|------------------|-----------|---------------------|------|
| Gateway          | IP1M12401 | Last battery status | 100% |
| Connected device |           |                     |      |
| Device Name      |           | Device Code         |      |
| BP5              |           | 004D320C8866        |      |
| 1000M            |           | C1BC779F0938        |      |

Fig 27

**Icon buttons** can also be provided (Fig 26 **A**). Once clicked, the device details information is displayed in a pop-up window (Fig 27). In the example provided all the wearable devices that were connected to the Gateway App during the last communication of the Gateway with the system are listed.



The information provided in this view can be not updated if the kit is disconnected from the network and at the same time one or more wearables are detached from the gateway.

## 10. User-Location Configuration




Fig 28

It is possible to associate users to one or more locations using the functionality **User-Location Configuration**. This procedure can be performed only by users with specific permissions. To do that,

- Click on three dots button (Fig 28 A) on Smart Monitor Web command bar and then on **Configure User-Location** button (Fig 28 B). The following window opens:




Fig 29

A grid is provided with the current associations between enabled locations and the active users. Locations can be shared between different users. The **"View by user"** (Fig 29 A) tab is selected by default and it shows, for every user, the list of associated locations (Fig 29 C). It is also possible to view the association by location, by selecting the **"View by location"** (Fig 29 D): for every location, the list of associated users is displayed. The existing associations can be removed:

- clicking the **Remove** button (Fig 29 B) provided on the right of each entry.
- While to configure a new association or to edit an existing one:

- click on **Create New** button placed above the grid.

The following window is displayed:

Fig 30

- Select a user on the "User" drop-down menu (Fig 30 A).
- Select the list of location(s) to be associated to the selected user on the "Location" drop-down menu (**Error! Reference source not found.B**).
- Click **Save** to confirm (Fig 30**Error! Reference source not found.C**).

From the version 8.1 of Digistat Suite package, the association between users and locations can be also configured through the Configurator Web > System Configuration > Locations .

## 11. Patient admission mode

| Location | Slot | Patient  | HR (bpm)   | RR (bpm)     | Skin Temp. (°C) | SpO2 (%) | NBPd (mm Hg) | NBPs (mm Hg) | Priority | INFO |
|----------|------|--|------------|--------------|-----------------|----------|--------------|--------------|----------|------|
| ICU      | 2    | Cartonio Tonio<br>T0002 Age: 33 Years<br>Sex: d  | 63<br>now  | 21.49<br>now | 89<br>now       | N/A      | 77<br>now    | 134<br>now   | N/A      | N/A  |
| ICU      | 3    | Cotogno Milo<br>M0003 Age: 23 Years<br>Sex: d    | 61<br>now  | 21.02<br>now | 97<br>now       | N/A      | 83<br>now    | 125<br>now   | N/A      | N/A  |
| ICU      | 4    | Lina Fata<br>F0004 Age: 34 Years<br>Sex: f       | 70<br>now  | 19.97<br>now | 24<br>now       | N/A      | 87<br>now    | 123<br>now   | N/A      | N/A  |
| ICU      | 5    | Lucente Amina<br>A0005 Age: 22 Years<br>Sex: f   | 63<br>now  | 21.71<br>now | 87<br>now       | N/A      | 87<br>now    | 115<br>now   | N/A      | N/A  |
| ICU      | 6    | Rosaropa Strega<br>S0006 Age: 21 Years<br>Sex: f | 71<br>now  | 23.07<br>now | 95<br>now       | N/A      | 87<br>now    | 131<br>now   | N/A      | N/A  |
| ICU      | 8    | Manno Orco<br>O0008 Age: 36 Years<br>Sex: d      | N/A<br>now | 23.76<br>now | N/A<br>now      | N/A      | N/A<br>now   | N/A<br>now   | N/A      | N/A  |
| ICU      | 1    | Fosco Lupo<br>L0001 Age: 40 Years<br>Sex: d      | 70<br>now  | 20.32<br>now | 87<br>now       | N/A      | 85<br>now    | 119<br>now   | N/A      | N/A  |

Fig 31

Smart Monitor Web can be configured to manage patients' admission and discharge instead of devices assignment and detachment.

In these cases, the devices assignment and detachment procedures are performed on an external system, depending on the specific choices of the healthcare organization.

When Smart Monitor Web is in "Patient admission" mode the Assign device/Detach device buttons on the Main Menu and on the toolbar drop down menu are replaced by the Admit Patient (Fig 31 **A**)/Discharge Patient (Fig 31 **B**) buttons.

## 11.1 Patient admission

To admit a patient:

- Click the **Admit Patient** button on the Smart Monitor Web command bar.

The following screen opens:

The screenshot shows the 'Admit Patient' window. It has a title bar 'Admit Patient' with a close button. The main area is divided into two sections. The left section is titled '1 SELECT PATIENT' and contains a text input field 'Insert Patient Code:' with 'PT011' entered, a 'SEARCH' button, and a section titled '2 SELECT LOCATION' with a 'Select Location:' dropdown set to 'ICU2' and a 'Select Bed:' dropdown with a list of beds (1, 2, 3, 4, 5) where bed '2' is selected. The right section is titled 'FOUND PATIENT' and displays patient information: 'TEST, PATIENT PT', 'Born on N/A in N/A', and 'Code: PT011'. There are four red callout markers: 'A' points to the 'Insert Patient Code:' field, 'B' points to the 'Select Bed:' dropdown, 'C' points to the 'FOUND PATIENT' section, and 'D' points to the 'Admit patient' button at the bottom right.

Fig 32

- Insert the **Patient Code** in the field indicated in **Error! Reference source not found. A** and click on **SEARCH** button.

If the patient is found, and they are not already admitted, the patient data is displayed on the right (**Error! Reference source not found. C**).

- Select the destination **Location** and **Bed** (**Error! Reference source not found. B**).

The screenshot shows a 'Patient admission confirmation' dialog box. It has a title bar 'Patient admission confirmation'. The main text asks 'Is all data correct? Do you really want to admit the patient?'. At the bottom, there are two buttons: 'CONFIRM' and 'CANCEL'.

Fig 33

- Click the **Admit patient** button on the bottom-right corner of the screen (**Error! Reference source not found. D**), then confirm the patient admission (Fig 33).

## 11.2 Patient Discharge

To discharge a patient:

- Click the **Discharge Patient** button on the Smart Monitor Web command bar (Fig 31 **B**).

The following screen opens:

Fig 34

- Insert the **Patient Code** in the field indicated in Fig 34 **A** and click on **SEARCH** button. If the patient is found the patient data is displayed on the right (Fig 34 **B**).

Fig 35

- Click on the **Discharge Patient** button (Fig 34 **C**) then confirm the patient dischargement (Fig 35).