multiAssist-module

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Versie 5.0

2024-03-25

Inhoud

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

The 'Move Patient' dialog box has a teal header. It contains two dropdown menus: 'Location *' with 'ICU' selected, and 'Bed *'. At the bottom are 'MOVE' and 'CANCEL' buttons.

Fig 14

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The 'NEW PATIENT' form has a teal header. It includes fields for 'Family Name *', 'Given Name *', 'Initials', 'Code *', 'Birth Date' (with a date picker), and 'Sex'. There is a 'Notes' text area. A red location pin icon with the letter 'A' is placed over the 'OPEN ADMISSION' button, which is highlighted with a red rectangle. At the bottom are 'CREATS' and 'CANCEL' buttons.

Fig 15

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1.3.1. Voorschriftformulier

Patient information: Weight 85 kg, Height 180 cm, BMI 26.23, Body Surface Area 2.06, Main Diagnosis, Medications, Allergies.

BGA: spO2 212 mmHg, pO2 100 mmHg, pH 98, cCa++ 48.

Choose anticoagulation: Heparin, CiCa.

Treatment mode: Post-CVVHDF, Pre-CVVHDF, CVVHD, Pre-CVVH, Post-CVVH, Pre-Post CVVH.

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Choose anticoagulation: Heparin, CiCa.

Treatment mode: Post-CVVHDF, Pre-CVVHDF, CVVHD, Pre-CVVH, Post-CVVH, Pre-Post CVVH.

Choose anticoagulation: Heparin, CiCa.

Treatment mode: Ci-Ca PostCVVHDF, Ci-Ca CVVHD.

Fig 30

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
Expected serum bicarbonate concentration with Ci-Ca postCVVHDF in relation to blood and dialysate flow: Simulation was performed with the following constant treatment parameters: Citrate dose 5.0mmol/L, Calcium dose 1.7 mmol/L, Net Ultrafiltration 100mL/h, Substitution flow = 17% of Blood flow (based on a numerical model - for orientation)

Fig 43

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- KDIGO recommendation: 20-25 ml/kg/h

Treatment machine dose set
34 ml/kg/h

- In order to achieve a delivered dose of 20-25 ml/kg/h, it is likely that the prescription will need to be in the range of 25-30 ml/kg/h

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36.5 ml/kg/h

- Uncommon prescription range

- In order to achieve a delivered dose of 20-25 mL/kg/h, it is likely that the prescription will need to be in the range of 25-30 mL/kg/h

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Set Disposables

multiAssist Prescription DASHBOARD EDIT DELETE ascom

Fig 46

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Pressures	Heparin	Anticoagulation	UF rate	Next operator action	History
				Filtrate bag change.	0 minute
				Dialysate bag change.	2 hours 20 minutes
				Citrate bag change.	6 hours 35 minutes _Change bag_
				Ca bag change	7 hours 45 minutes _Change bag_

Fig 75

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1.4.7. Tabblad "Geschiedenis"

Pressures	Heparin	Anticoagulation	Uf rate	_Next action_	History
11/20/23, 2:40 PM	Created				ADMIN
11/20/23, 2:43 PM	PostFilterCalciumEdited			Post-filter calcium set as 0.12 mmol/L at 11/20/2023 13:45:00.	ADMIN
11/20/23, 2:45 PM	PostFilterCalciumEdited			Post-filter calcium set as 0.06 mmol/L at 11/20/2023 13:39:00.	ADMIN
11/20/23, 2:54 PM	PostFilterCalciumEdited			Post-filter calcium set as 0.31 mmol/L at 11/20/2023 13:49:00.	ADMIN
11/20/23, 2:54 PM	PostFilterCalciumEdited			Post-filter calcium set as 0.46 mmol/L at 11/20/2023 13:51:00.	ADMIN
11/20/23, 3:37 PM	Information			Level in bubble catcher has dropped	
11/20/23, 3:41 PM	Information			Low access pressure alarm	
11/20/23, 3:42 PM	Information			DiaF/BF ratio < 25 %	
11/20/23, 3:43 PM	Information			DiaF/BF ratio < 25 %	

Fig 80

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Prescription

Prescribed

Missing prescription click on button above to create a prescription with the data from device

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The screenshot displays the multiAssist Module interface, which is organized into several sections:

- Patient information:** A sidebar on the left containing patient details such as Weight (78 kg), Height (176 cm), Bmi (25.18), Body surface area (1.95), Medications (No), Allergies (No), and various blood gas and electrolyte values (pH, pO2, sO2, fO2Hb, fCOHb, fHHb, fMetHb, cCl-, ctO2, ctO2, Acid base Status).
- Choose anticoagulation:** A section with buttons for Heparin and CiCa.
- Treatment mode:** A section with buttons for Ci-Ca postCVVHDF and Ci-Ca CVVHD.
- Dialysis parameters:** A central area with a 'RESET DEFAULT' button and various adjustable parameters:
 - DialysateFlow: 2000 ml/h
 - BloodFlow: 110 ml/min
 - UltrafiltrationRate: 10 ml/h
 - SubstituteFlow: 600 ml/h
 - Temperature: 38 °C
- Renal Dose:** A section on the right showing:
 - Default renal dose: 51 ml/kg/h
 - Target renal dose: 43 ml/kg/h
 - Predicted renal dose: 38.7 ml/kg/h
- Anticoagulation settings:** A section at the bottom with buttons for Heparin and Uf goal, and fields for CitrateDose (5 mmol/l) and CalciumDose (1.7 mmol/l).

Fig 87

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1.6. Appendix

This document provides the references for data and formulas used by multiAssist module of Digistat Care product. The formulas have been validated by Fresenius team as detailed in the [“CCTS Formulas 20231207 ACDM-585-14838”](#) document.

1.6.1. Definitions

1.6.1.1. Parameter Definitions

Definitions of parameter applied in formulas to calculate the Treatment machine dose set. All parameters are those set by the user except where otherwise specified.

Term	Unit	Description
RD	ml/h/kg	Renal dose
UF	ml/h	Ultrafiltration rate
HF	ml/h	Heparin flow
BF	ml/min	Blood flow
DF	ml/h	Dialysate flow
preSF	ml/h	Pre-substitute flow
postSF	ml/h	Post-substitute flow
Ci _D	mmol/l	Citrate dose
Ca _D	mmol/l	Calcium dose
CiF	ml/h	Citrate Flow ⁴
CaF	ml/h	Calcium flow ⁵
Ci _C	mmol/l	Citrate concentration solution ¹
Ca _C	mmol/l	Calcium concentration solution ²
CaS _C	mmol/l	Calcium substitute concentration solution ³
W	kg	Weight of the patient

¹ Ci_C from system-option “CitrateConcentration”

² Ca_C from system-option “CalciumConcentration”

³ CaS_C from system-option “CalciumSubstituteConcentration”

⁴ CiF Calculated

⁵ CaF Calculated

1.6.1.2. Machine parameters

List of parameters acquired by MultiFiltratePro. For each parameter:

- Name: name of the parameter when used into formulas
- ParID: Digistat StandardParameterID
- UomID: Digistat StandardUnitofMeasureID
- Description: a textual description of the parameter
- FMC Code: code of the parameter in the Fresenius MultifiltratePro protocol.

Name	ParID	UomID	Description	FMC Code
BF	5500	1306	Blood Flow set	QFOYSY
	5110	1306	Actual Blood Flow	QPOYSY
St	5114	99999	Treatment status code	GFZSON
	5113	36	Filter lifetime	QPZSSM
	5087	36	Treatment time, running total	QPZSGG
Ci-Ca	5532	99999	Ci-Ca Active	QFPPNP
	5104	99999	Treatment Start	GFGFGE
C	5070	99999	Therapy Type	GFZSGU
	5518	1307	Substitution rate set	QFZSNE
DF	5503	1307	Dialysate Flow set	QFZSQE
UF	5517	1307	UF rate set	QFZSRR
	5072	1307	Citrate Flow	QPZPVS
CaF	5073	1307	Calcium Flow	QPZPNS
UFR	5020	1002	Ultrafiltration Ratio	QPHSOS
UV	5079	1302	Effective total withdrawal total	QPZSRG
SV	5083	1302	Substitute-Total total	QPZSFG
DV	5077	1302	DF-Total balance total	QPZSQG
CiV	5095	1302	Citrate volume, total	QPZSVT
CaV	5093	1302	Calcium volume, total	QPZSNT
HV	5089	1302	Anticoagulation, continuous total	QPZSXG
	5508	1307	Heparin Continuous Rate set	QFURCE
	5516	22	Dialysate temperature set	QFQVGC
	5520	1642	Citrate Dose set	QFZPVQ
	5521	1642	Calcium Dose set	QFZPNQ
	5096	3	Weight of scale 1	QPZSJN
	5097	3	Weight of scale 2	QPZSJO
	5098	3	Weight of scale 3	QPZSJP
	5099	3	Weight of scale 4	QPZSJQ
	7583	99999	Device Serial number	GFFEAO
SF	5518	1307	Substitution rate setting	QFZSNE
	5107	99999	Treatment end	GFGRAQ
	5108	99999	Treatment Pause Begin	GFGCNO
Sc	5109	99999	Treatment Pause End	GFGCNR
	5533	99999	Device Status	GFZSFP
	5001	1401	Access Pressure	QPNECE
	5004	1401	Return Pressure	QPIRCE
	5006	1401	Transmembrane Pressure	QPGZCI
	5534	99999	Alarms and warnings on the device	QPZSRE
	5112	1401	Pre-Hemofilter pressure	QPZSCU
	5003	1401	Effluent Pressure	QPZSCS

1.6.2. Formulas

1.6.2.1. DF/BF ratio

Dialysate/Blood flow ratio is defined as:

$$r = \frac{DF}{BF * 60}$$

Where:

$$r = \frac{UF + postSF + CaF}{BF * 60}$$

BF: Blood flow as set by the user

DF: Dialysate flow as set by the user

1.6.2.2. UF/BF ratio

Ultrafiltration/Blood flow ratio is defined as:

Where:

Parameter	Unit	Description
r		ratio
UF	ml/h	Ultrafiltration rate as set by the user
BF	ml/min	Blood flow as set by the user
postSF	ml/h	Post-substitute flow as set by the user (equals to 0 if the therapy doesn't have the "Post-substitute flow" parameter)
CaF	ml/h	Calcium flow as calculated (equals to 0 if therapy is not a Ci-Ca therapy)

1.6.2.3. Calcium Flow

Calcium Flow is defined as:

$$CaF = \frac{Cad * (DF + UF + CiF + HF + postSF) - CaSc * postSF}{Cac - Cad}$$

1.6.2.4. Citrate Flow

Citrate Flow is defined as:

$$CiFp = \frac{Cid * BF * 60}{Cic}$$

1.6.2.5. PreDilutionFactor Prescription

PreDilutionFactor used in prescription is defined as:

$$Pfp = \frac{(BF * 60)}{(BF * 60) + CiF + HF + preSF}$$

1.6.2.6. Machine Renal Dose Set

Renal dose calculation depends on the therapy as follows:

- For a **Post-CVVHDF** therapy, treatment machine dose set is calculated by this formula:

$$RD = \frac{UF + HF + DF + postSF}{W}$$

- For a **Pre-CVVHDF** therapy, treatment machine dose set is calculated by this formula:

$$RD = Pfp * \frac{UF + HF + DF + preSF}{W}$$

- For a **CVVHD** therapy, treatment machine dose set is calculated by this formula:

$$RD = \frac{UF + HF + DF}{W}$$

- For a **Pre-CVVH** therapy, treatment machine dose set is calculated by this formula:

$$RD = Pfp * \frac{UF + HF + preSF}{W}$$

- For a **Post-CVVH** therapy, treatment machine dose set is calculated by this formula:

$$RD = \frac{UF + HF + postSF}{W}$$

- For a **Pre-Post-CVVH** therapy, treatment machine dose set is calculated by this formula:

$$RD = Pfp * \frac{UF + HF + preSF + postSF}{W}$$

- For a **Ci-Ca postCVVHDF** therapy, treatment machine dose set is calculated by this formula:

$$RD = Pfp * \frac{UF + HF + DF + postSF + CiF + CaF}{W}$$

- For a **Ci-Ca CVVHD** therapy, treatment machine dose set is calculated by this formula:

$$RD = Pfp * \frac{UF + HF + DF + CiF + CaF}{W}$$

1.6.2.7. Substitute Flow (SFc)

Substitute Flow (SFc) is needed to determine the current therapy ongoing on the MultiFiltratePro.

It is calculated as:

$$SF_c = \frac{UF_R * BF * 60}{100} - UF - CF$$

1.6.2.8. Treatment Time (TT)

The duration of treatment [TT] shall be calculated after determining the start date of treatment. The start date corresponds to the start date of the blood pump, indicated by the duration of use of the filter (parameter **FLT** in minutes). This parameter is only updated every 10 minutes, so the system needs to detect the change of the value to calculate the start date. The interception must be in the interval of 1 minute.

Example 1:

Dataset time (date part is not relevant)	Filter lifetime
10:33	0
10:34	0
...	...
10:44	0
10:45	10
10:46	10
...	...

←Change detected.

The start time of the blood pump is considered as: $10:45 - 10\text{min} = 10:35$

Example 2:

Dataset time (date part is not relevant)	Filter lifetime
10:33	0
10:34	0
...	...
10:45	0
10:47	10
...	...
10:54	10
10:55	20
10:56	20
...	...

←the interval is 2 minutes: no interception


←Change detected

The start time of the blood pump is: $10:55 - 20\text{min} = 10:35$

1.6.2.9. Therapy Code (Tc)

Code of the therapy running on the device is retrieved by using the following diagram. It requires the following calculated value first:

- Substitute Flow calculation (SFc)

1.6.2.10. Predilution Factor (Pf)

Predilution factor is used by "PRE" therapies and it's calculated as follows:

- This formula is used for all "PRE-only" therapies and Ci-Ca (Pre-CVVHDF, Pre-CVVH, Ci-Ca postCVVHDF, Ci-Ca CVVHD):

$$\frac{(BF * 60)}{(BF * 60) + CiF + HF + SF}$$

- This formula is used for "PRE-POST" therapies (Pre-Post CVVH):

$$\frac{(BF * 60)}{(BF * 60) + CiF + HF + DF}$$

1.6.2.11. Effluent Rate (ER)

The Effluent Rate is calculated only when all following occur:

- the device status is "Treatment" [**St** = "Treatment"]
- the device balance status is not stopped [**Sc** = 00 or **Sc** == 01]
- the effective treatment code is determined [**Tc**]
- Treatment time is determined [**Tt**]

This formula is used for the following therapies (all "PRE" therapies and CiCa):

- Pre-CVVHDF, Pre-CVVH, Pre-Post CVVH, Ci-Ca postCVVHDF, Ci-Ca CVVHD:

$$\frac{CiV + CaV + SV + DV + UV + HV}{\frac{TT}{60}} * pf$$

- Other therapies (Post-CVVHDF, CVVHD, Post-CVVH) adopt the following:

$$\frac{CiV + CaV + SV + DV + UV + HV}{\frac{TT}{60}}$$

1.6.2.12. Average Renal Dose

Average Renal Dose [RD] calculation requires the Effluent Rate [ER] and the patient weight [W]. It's calculated as follows:

$$RD = \frac{ER}{W}$$

1.6.2.13. Bags Calculations – Filtrate bag full

Calculations to Filtrate bag full is performed as follow (Note: all calculations consider a density of 1 g/ml)

Calculation can be performed only if all the following are satisfied:

- There is a prescription saved;
- Prescription has been saved with a disposable of type "Filtrate";
- Selected disposable has a volume set;
- There is an ongoing therapy on the machine;
- The machine is in status "Treatment";
- The balancing is active (see "4.8 Therapy Code (Tc)")

Filtrate Bag Target (Fbt) is retrieved by "volume" property of selected disposable.

FbT must take into account the Tara of the bag. Tara is retrieved as described in 1.6.3.3 Bags Volume/Tara relationship.

Filtrate bag Target effective (FbTe) = FbT – Tara

Calculation of Filtrate bag Speed (FbS):

Speed is calculated taking two consecutive dataset of parameters QPZSJP + QPZSJQ (weights on scales 3 and 4)

FbW (Filtrate bag weight) = QPZSJP + QPZSJQ

Given the following:

DateTimeClinical	FbW
T1 (last before current)	FbW1
T2 (current)	FbW2

FbS (Filtrate bag speed) is calculated as:

$$FbS = \frac{FbW2 - FbW1}{T2 - T1 (*)}$$

(*) = in minutes.

FbS must be greater than 0. If less or equal to 0, calculation cannot be performed.

Remaining time is calculated as:

$$Time = \frac{FbTe - FbW2}{FbS}$$

If resulting time is negative, it shall be considered as 0.

Resulting time is rounded with steps of 5 minutes (example: 13 minutes will become 10 minutes).

1.6.2.14. Bags Calculations – Dialysate/Substitute bags empty

Calculations to Dialysate/Substitute bag empty is performed as follow (Note: all calculations consider a density of 1 g/ml)

Calculation is performed only when all the following are met:

- There is an ongoing therapy on the machine;
- The machine status is "Treatment";
- Balancing is active (see 1.6.2.9 Therapy Code (Tc)).

There could be different calculation depending on the therapy that is running on the device retrieved as described in 1.6.2.9 Therapy Code (Tc).

For all therapies except CiCA Post CVVHDF calculation is performed on the sum of the scales 1 and 2 given by parameters QPZSJN and QPZSJO (see 1.6.3.4 Scale-Bags-Therapy relationship, as reference)

DialSub bag weight (DSbW) = QPZSJN + QPZSJO

Target weight is considered always as 900 grams (we assume that the max load of 20 kgs will apply on balances 1 and 2)

DialSub bag Target (DSbT) = 900 g

Speed is calculated taking two consecutive datasets as follows:

DateTimeClinical	DSbW
T1 (last before current)	DSbW1
T2 (current)	DSbW2

DSbS (DialSub bag speed) is calculated as:

$$DSbS = \frac{DSbW2 - DSbW1}{T2 - T1 (*)}$$

(*) = in minutes.

Due to the nature of the calculation, DSbS is usually negative. If DSbS is equal or greater than 0, time calculation will not be performed.

Remaining time is calculated as:

$$Time = \frac{DSbW2 - DSbT}{-DSbS}$$

Resulting time is rounded with steps of 5 minutes (example: 13 minutes will become 10 minutes)

Output message shall be produced according to the running therapy type as described in 1.6.3.4 Scale-Bags-Therapy relationship. As an example, if a CVVHD therapy is running the message shall refer to the emptying of Dialysate.

For CiCA Post CVVHDF therapy calculation is performed on each scale given by parameters QPZSJN and QPZSJO (see 1.6.3.4 Scale-Bags-Therapy relationship, as reference)

Disalysate bag weight (DbW) = QPZSJN
Substitute bag weight (SbW) = QPZSJO

Target weight is considered always as 500 grams (we assume that the max load of 10 kgs on scale 1 and another 10 Kgs on scale 2 will apply)

For Dialysate bag the following applies (for Substitute the calculation is exactly the same except the data source, SbW = QPZSJO, and the returned message shall refer to substitute)

Dialysate bag Target (DbT) = 500 g

Speed is calculated taking two consecutive datasets as follows:

DateTimeClinical	DbW
T1 (last before current)	DbW1
T2 (current)	DbW2

DbS (DialSub bag speed) is calculated as:

$$DbS = \frac{DbW2 - DbW1}{T2 - T1 (*)}$$

(*) = in minutes.

Due to the nature of the calculation, DbS is usually negative. If DbS is equal or greater than 0, time calculation will not be performed.

$$Time = \frac{DbW2 - DbT}{- DbS}$$

If resulting time is negative, it shall be considered as 0.

Resulting time is rounded with steps of 5 minutes (example: 13 minutes will become 10 minutes).

1.6.2.15. Bags Calculations – Citrate/Calcium bag empty

Calculations to Citrate and Calcium bags empty is performed as follows (Note: all calculations consider a density of 1 g/ml)

Calculation can be performed only if all the following are satisfied:

- There is a prescription saved;
- Prescription has been saved with a disposable of type "Citrate"/"Calcium";
- Selected disposable has a volume set;
- There is an ongoing therapy on the machine;
- The machine is in status "Treatment";
- The balancing is active (see 1.6.2.9 Therapy Code (Tc)).

Calculation relies on:

- Therapy treatment time (see 1.6.2.8 Treatment Time (TT));
- Manual action on multiAssist to record bag change;
- Citrate flow (QPZPVS).

Calculus

BagFlow = QPZPVS for Citrate Flow (QPZPNS for Calcium Flow) retrieved from the machine

BagVol = the volume of the Citrate bag (or Calcium Bag) as retrieved by disposable selected in prescription.

ChangeT = the most recent date/time from:

- Therapy start time (ThStart), calculated as:
ThStart = Current date/time – Treatment Time (see 1.6.2.8 Treatment Time (TT));
- A manual action of “bag change” has been recorded in multiAssist.

Remaining time for Citrate (Calcium) bag empty is calculated as:

Given DeltaT as:

$$DeltaT = (DateTimeCurrent - ChangeT)$$

BagFlow must be greater than 0. If 0 or negative, Time calculation will not be performed.

Remaining time is calculated as:

$$Time = \frac{(BagVol - BagFlow * DeltaT)}{BagFlow}$$

If resulting time is negative, it shall be considered as 0. Resulting time is rounded with steps of 5 minutes (example: 13 minutes will become 10 minutes).

1.6.3. Tables

1.6.3.1. Prescription parameters limits and default values

The following table defines the parameter values for each therapy: the min-max values and the default value (in brackets). All the parameters with default value are required, except the heparin flow.

	Blood flow in ml/min	Temperature in °C	Dialysate flow in ml/h	Pre-substitute flow in ml/h	Post-substitute flow in ml/h	Ultrafiltration rate in ml/h	Heparin flow in ml/h	Citrate dose in mmol/l	Calcium dose in mmol/l
Resolution →	10	0,5	10	10	10	10	0,1	0,1	0,1
Pre-CVVHDF	10-500 (100)	35-39 (38)	600-4800 (1000)	600-4800 (1000)		0-990 (0)	0,5-25,0 (off:0 – on: 0,5)		
Post-CVVHDF	10-500 (100)	35-39 (38)	600-4800 (1000)		600-4800 (1000)	0-990 (0)	0,5-25,0 (0,5)		
CVVHD	10-500 (100)	35-39 (38)	600-4800 (2000)			0-990 (0)	0,5-25,0 (0,5)		
Pre-CVVH	10-500 (100)	35-39 (38)		600-4800 (1000)		0-990 (0)	0,5-25,0 (0,5)		
Post-CVVH	10-500 (100)	35-39 (38)			600-4800 (1000)	0-990 (0)	0,5-25,0 (0,5)		
Pre-post CVVH	10-500 (100)	35-39 (38)		600-4800 (1000)	600-4800 (1000)	0-990 (0)	0,5-25,0 (0,5)		
Ci-Ca postCVVHDF	10-200 (100)	35-39 (38)	600-4800 (2000)		600-2400 (1000)	0-990 (0)	0,5-25,0 (0,5)	2,0-6,0 (4,0)	0,1-3,0 (1,7)
Ci-Ca CVVHD	10-200 (100)	35-39 (38)	600-4800 (2000)			0-990 (0)	0,5-25,0 (0,5)	2,0-6,0 (4,0)	0,1-3,0 (1,7)

1.6.3.2. Prescription parameters limits and default values

For a running therapy, user can switch prescription according to the following table:

Initial therapy defined	Therapies allowed
Pre-CVVHDF	Post-CVVHDF, CVVHD, Pre-CVVH, Post-CVVH
Post-CVVHDF	Pre-CVVHDF, CVVHD, Pre-CVVH, Post-CVVH
CVVHD	Pre-CVVHDF, Post-CVVHDF, Pre-CVVH, Post-CVVH
Pre-CVVH	Pre-CVVHDF, Post-CVVHDF, Post-CVVH, CVVHD
Post-CVVH	Pre-CVVHDF, Post-CVVHDF, Pre-CVVH, CVVHD
Pre-post CVVH	Pre-CVVH, Post-CVVH
Ci-Ca postCVVHDF	Pre-CVVHDF, Post-CVVHDF
Ci-Ca CVVHD	CVVHD

1.6.3.3. Bags Volume/Tara relationship

The following table represents the relationship between the applied weight on balances and the tara:

Applied weight in g	Tara in g
0 < weight <= 7000	300
7000 < weight <= 12000	500
12000 < weight <= 17000	700
weight <=17000	900

Weight-Volume relationship is considered by applying a default density of 1 g/ml.

1.6.3.4. Scale-Bags-Therapy relationship

The following table shows the relationship between device scales, bag types and therapies.

Bag	Therapy							
	Pre-CVVHDF	Post-CVVHDF	CVVHD	Pre-CVVH	Post-CVVH	Per-Post-CVVH	CiCa post-CVVHDF	CiCa CVVHD
Dialyate			Scale 1 and 2				Scale 1	Scale 1 and 2
Substitute				Scale 1 and 2	Scale 1 and 2	Scale 1 and 2	Scale 2	
Dial/Sub	Scale 1 and 2	Scale 1 and 2						
Ca							Left pole	Left pole
Citrate							Left pole	Left pole
Filtrate	Scales 3 and 4							