



## Switch Modules/Pull Cord Switched on for security and service

Plugged in and switched on. For whatever reason, when they hit the button or pull the cord, your patients will know that their message will be communicated in a split second and that all necessary assistance will soon be on its way.

Specially developed to enhance any hospital or nursing home environment, the ascom teleCARE M Switch Modules/Pull Cord provide your patients with security and service, and you with versatility, reliability and value for money.

### **Value through versatility...**

Freely combined in any teleCARE M system, up to three switches (each containing one, two or three function-button combinations) are installed in the patient's room or bathroom to support each patient's individual requirements. Everything from normal call, assistance call, nurse presence, cancel and bed-light control can be catered for and simply modified as and when required.

### **...reliability and smart design**

The reliability of the system is second to none and downtimes can be kept to an absolute minimum. For example, a Safe Release Socket releases the plug if excessive strain is applied to the cable connecting the teleCARE M patient handset to the switch module. This reduces the number of damaged plugs and broken cables, and significantly saves repair costs and increases system reliability.

Further smart features such as LED indicated function buttons, integrated pull-cord switches and an ergonomic design, ensure that the ascom teleCARE M Switch Modules/Pull Cord are easy to use and maintain, and that patients are comfortable using them at all times.

## Features

### **Integrated pull-cord switches**

Available for bathroom and toilet use.

### **Safe Release socket**

Releases the plug and saves repairs.

### **LED indicated function buttons**

Dull light for night location, bright light during day.

### **Full flexibility**

Supports all combinations on bed level.

### **Single, double and triple backplates**

Easy installation and maintenance. Can be mounted over various international back boxes.

### **Ergonomic design**

Easy to clean and maintain.

## Basic Data

Size: 91 x 91 x 17 mm

Colour: White or light grey

## Accessories

Wet Area Cover for single or pull-cord switches

Surface Mounting Spacer

Pull cord – pull force to generate an alarm 2N

