

ascom teleCARE
Wireless Switch Modules



Up to three functions and a connection socket in a single wireless switch module

A unique combination of design and functionality characterises the teleCARE wireless switch modules. These extremely compact modules require no external power supply and can be mounted virtually anywhere. The teleCARE wireless switch modules have three function buttons and can even incorporate the unique Safe Release socket in a single module or a pull-cord call switch.

Each function button contains an LED which can be programmed to emit a bright light for up to 10 seconds after the circuit is activated. In addition to this visual reassurance signal, the teleCARE wireless switch modules also incorporate an optional reassurance buzzer. The compact design of the teleCARE wireless

switch modules and the specially designed mounting box means that they can easily and rapidly be installed in any suitable location. The teleCARE wireless switch modules can also be used in combination with the more conventional "hard-wired" teleCARE switch modules.

The unique teleCARE Safe Release socket is so compact that it can be integrated in the standard sized wireless switch modules. The Safe Release socket is for connection of the teleCARE patient handsets. Its revolutionary design eliminates the problems of damaged plugs and broken cables by releasing the plug if excessive strain is applied to the handset cable.

Main Features

- Attractive styling and an ergonomic design with a neutral grey colour which blends with any background.
- Available with three function buttons which are individually coloured to indicate their functions and each button has an LED.
- The unique teleCARE Safe Release socket can be incorporated in the standard size wireless switch module.
- The standard sized wireless switch module can incorporate a pull-cord call switch.
- Compact design and no need to penetrate into the wall means that the teleCARE wireless switch modules can be easily mounted almost anywhere.
- Fully compatible with the teleCARE hardwired switch modules.
- The wireless switch module contains a replaceable 9V battery which is permanently monitored and an automatic low-power alarm signal is transmitted when the battery output falls below the preset safe level.



3-button module with socket



3-button module



3-button module with pull-cord



Safe Release plug

ascom

Ascom Tateco AB
Wireless Solutions

P.O. Box 8783 SE-402 76 GÖTEBORG
Tel +46 31 55 93 00 Fax +46 31 55 20 31
www.ascom.com/ws

Specifications are subject to change without notice.