



## Ascom Telligence room devices — for bedside care and workflow efficiency

*Ascom Telligence room devices are peripheral devices and modules that help enhance staff efficiency and support patient safety measures. They include doorside room displays, bedside and doorside modules, an external input module for medical and non-medical devices, pull cords, a standalone speech module, and an input/output interface unit.*

### Key benefits

- Improve patient and staff satisfaction with customized bedside modules. Each button's color, text and icon can be changed to suit a wide range of needs.
- Enhance staff communication and collaboration. Doorside room displays show text messages that are legible from virtually anywhere in a room. Incoming voice calls are automatically prioritized. Speech modules let caregivers communicate with colleagues while remaining with the patient.
- All Ascom Telligence bedside devices now feature Ascom SafeConnect™ magnetic receptacles and plugs, helping to minimize equipment damage should the connector be accidentally disconnected, and sends an alert if there is a problem.
- Ensure uninterrupted care. Ascom Telligence room devices are 'hot swappable'.
- Help keep patients safe. Ascom Telligence room devices include a speech module that lets staff communicate from toilet areas without leaving the patient.
- Cost efficient. Modular concept lets you build tailor-made systems for specific rooms and requirements. Systems can be modified and expanded as needs and care models change.

Contact your nearest Ascom representative at [ascom.com](http://ascom.com) to learn how the Ascom Telligence Patient Response System can help improve efficiency and patient care.



## Key features, Ascom Telligence ROOM devices:

- RJ-45 connectors enable the use of pre-terminated cables, reducing the installation time compared to traditional connectors and proprietary cables.
- SafeConnect magnetic plugs. Absence of pins prevents damage to devices, even when plugs are forcefully removed.
- All Ascom Telligence room devices are constructed from durable materials that let them withstand repeated and thorough cleaning.
- Speech modules deliver crisp, high-quality voice quality.

### DEVICE LIST

	Customizable Button EU, NUDM3-HE		Speech Module, EU NUSPM-HE
	Customizable Button, Passive EU, NUDM3P-HE		Emergency Pull, UK NUEPS-HK
	Customizable Bedside Module EU, NUBM3-HE		Cardiac Pull, UK NUCPS-HK
	Bedside Extension Module NUBM3X-HE		IO Module, Active NUIOM-HG
	Customizable Bedside Module, Passive, EU NUBMP-HE		Pull cord IP44, EU NUPC-HE
	External Input Module EU, NUM2-HE		
	Room Display NURD-HE		Pull Cord White 3b Passive NIPC-W3P
	Room Display-Passive RMD		

*These room devices are used as components in Ascom Telligence solutions. Contact Ascom today to learn how Ascom Telligence can help hospitals and patients.*

**Ascom Holding AG**  
Zugerstrasse 32,  
CH-6340 Baar,  
Switzerland  
info@ascom.com  
Phone: +41 41 544 78 00  
ascom.com

**ascom**

The availability, configuration and technical specifications of Ascom products, services and solutions may vary from country to country. Please ask your nearest Ascom representative for further details.