

[ PRODUCT SHEET ]

**Product:** Ascom teleCARE patient handsets

From simple to sophisticated functionality

Ergonomic

Intuitive use

Contemporary design



## ASCOM teleCARE PATIENT HANDSETS – SUPERIOR SERVICE AND FUNCTIONALITY FOR CURE AND CARE ENVIRONMENTS



Wherever Ascom teleCARE patient handsets are installed, improved patient comfort, safety and peace of mind are always within reach.

Ascom teleCARE patient handsets are renowned for their ease of use and robustness. For 2012, they receive significant design and functionality improvements.

The new product line features four ergonomic handsets that offer patients an extensive range of features and capabilities, including: alarm and service calls, bed and room lighting, integrated flashlight, window shade/radio/TV controls, headset connection and nurse call speech.

### Functionality and design

All handsets feature color-coded and clearly marked buttons, with international symbols that identify their functionality.

The dominant feature on all handsets is the Call Button – a large, red button with a nurse symbol. Like other functions on the handset, the call button is permanently backlit to help patients locate it in the dark. When pressed, the light brightens to reassure patients that calls have been registered.

### Durable quality

All handsets are lightweight, compact, waterproof and can be dip-sterilized. The Safe Release plug, which connects each handset to the switch module, automatically releases if excessive strain is applied to the handset cable.

### RFID tag - an industry first!

To meet the demands of an ever-increasing inventory of products in a healthcare environment, Ascom teleCARE patient handsets now feature an integrated RFID tag. This tag can be used to simplify handset inventorying and sterilization routines.

[ PRODUCT SHEET: TELECARE PATIENT HANDSETS ]

| Ascom teleCARE patient handset comparison |         |                        |   |                          |  |
|---|---------|------------------------|---|--------------------------|--|
| Feature                                   | Variant | NIPH2-A1A<br>(Low end) | NIPH2-A3A (2,5m)<br>NIPH2-A3L (6m)<br>(Low end) | NIPH3-A7A<br>(Mid range) | NIPH3-AES<br>(High end)                    |
| Call button                               |         | ✓                      | ✓   | ✓                        | ✓  |
| Call reassurance led                      |         | ✓                      | ✓   | ✓                        | ✓  |
| Night orientation call button             |         | ✓                      | ✓   | ✓                        | ✓  |
| Service button                            |         | -                      | -   | ✓                        | ✓  |
| Service button reassurance led            |         | -                      | -   | ✓                        | ✓  |
| Night orientation service button          |         | -                      | -   | ✓                        | ✓  |
| Bed light control                         |         | -                      | ✓ (2 sources)                                   | ✓ (2 sources)            | ✓ (2 sources)                              |
| Orientation bed light buttons             |         | -                      | ✓   | ✓                        | ✓  |
| Sunblind control                          |         | -                      | -   | ✓ (2 sources)            | ✓ (2 sources)                              |
| Flashlight                                |         | -                      | -   | ✓                        | ✓  |
| Radio and TV sound                        |         | -                      | -   | -                        | Stereo on headset,<br>Mono handset speaker |
| Radio control                             |         | -                      | -   | -                        | FM radio                                   |
| Number radio channels                     |         | -                      | -   | -                        | 8 presets, or Scanned FM stations          |
| TV control                                |         | -                      | -   | -                        | ✓ infrared                                 |
| Support for IR code addressable TV sets   |         | -                      | -   | -                        | ✓  |
| Number of IR protocols                    |         | -                      | -   | -                        | 8  |
| Number of TV channels                     |         | -                      | -   | -                        | Depending TV set                           |
| Volume up/down                            |         | -                      | -   | -                        | ✓  |
| Channel up/down                           |         | -                      | -   | -                        | ✓  |
| Headset support                           |         | -                      | -   | -                        | ✓ (stereo), 3.5 mm jack                    |
| Nurse call speech standard                |         | -                      | -   | -                        | ✓ via internal speaker                     |
| Nurse call speech optional                |         | -                      | -   | -                        | Via bedside speaker                        |
| Telephone function                        |         | -                      | -   | -                        | -  |
| RFID tag in handset                       |         | -                      | -   | ✓                        | ✓  |
| Parameter configuration via WIN-PDM       |         | -                      | -   | ✓                        | ✓  |
| Bi-color design                           |         | ✓                      | ✓   | ✓                        | ✓  |
| Ingress protection                        |         | IP65                   | IP65  | IP65                     | IP65                                       |
| Safe release connector                    |         | ✓                      | ✓   | ✓                        | ✓  |
| Anti slip strip at the backside           |         | ✓                      | ✓   | ✓                        | ✓  |
| Dimensions (WxHxD) in mm                  |         | 54x151x20              | 54x151x20                                       | 54x151x20                | 54x151x20                                  |
| Cable length                              |         | 2,5 Meter              | 2,5 and 6 Meter                                 | 2,5 Meter                | 2,5 Meter                                  |
| <b>System related functions</b>           |         |                        |   |                          |  |
| teleCARE M and teleCARE IP compatible     |         | ✓                      | ✓   | ✓                        | ✓  |
| Red button with own call ID               |         | ✓ <sup>1</sup>         | ✓ <sup>1</sup>                                  | ✓ <sup>1</sup>           | ✓ <sup>1</sup>                             |
| Alarm when plugged out with own ID        |         | ✓ <sup>1</sup>         | ✓ <sup>1</sup>                                  | ✓ <sup>1</sup>           | ✓ <sup>1</sup>                             |
| Set-up parameters for project tailoring   |         | -                      | -   | ✓                        | ✓  |
| Firmware version in system diagnostics    |         | -                      | -   | ✓                        | ✓  |
| <b>Accessories</b>                        |         |                        |   |                          |  |
| Parking bracket                           |         | ✓                      | ✓   | ✓                        | ✓  |
| Linen clip                                |         | ✓                      | ✓   | ✓                        | ✓  |
| Cable replacement possible                |         | -                      | -   | ✓                        | ✓  |
| TV audio interface                        |         | -                      | -   | -                        | ✓  |

1 teleCARE IP only