



# Ascom teleCARE IP Nurse Call Waterproof Call Points

## Introducing the NIWP Wet Area Nurse Call Point

The NIWP call points are IP65 and AS3811 compliant modules with wet area suitable housings, with a great ergonomic, durable design. The call point has been designed to utilise the existing teleCARE IP Doorside Module (NIDM), and it is mounted to the wall via the standard Clipsal C clip. The NIWP are available in three variants. The function of the buttons can be changed to different configurations.

### Key Features

- LED light for each button to indicate the activated condition and night time guiding
- Waterproof device (IP65)
- Designed to be used in wet areas
- Simple installation to surface mount or to the wall box
- Programmable buttons

### Available Options

The NIWP are available in the following variants:

Option 1. Two buttons, nurse call, and cancel/nurse presence.

Option 2. Three buttons, nurse call, staff call, and cancel/nurse presence.

Option 3 - Three Buttons, emergency call, staff call and cancel/nurse presence.



### Ascom

41-43 Bourke Rd  
Alexandria, NSW  
Australia

Phone: 02 9698 9000  
[www.ascom.com/au](http://www.ascom.com/au)