




Ascom DURAlarm Product Sheet

Designed to promote alarm accuracy and enable rapid response to staff duress alarm, DURAlarm can receive alarms and messages from any source and deliver that information as a text message to any destination.



DURAlarm is the alarm module of Ascom DURAsuite Solution, which enables LAN and/or WiFi-connected users with a web browser to view in real-time, via text and graphical information, all alarms from wireless duress handhells, wired duress buttons and management systems.

When acknowledging an alarm, the user can add case notes for management purposes. All alarms can be sent to a wide variety of devices simultaneously so that the emergency response team can be instantly notified and can effectively manage the crisis situation.

All alarms are time-stamped and logged for reporting purposes.

Ascom DURAlarm Features:

- Real-time alarm management
- Map presentation with alarm location from the basestation, Locator and fixed device locations
- Map or Text alarm pop-ups
- Accept / Reject tagging and alarm escalation
- Text and speech alarm notifications on handheld devices
- Unlimited devices and device groups
- Unlimited users and user zones
- Map watermarking for easy map recognition
- Resilient LINUX based appliance with no moving parts
- No software installations required - fully web-based
- Logging of all alarm activity, including events and audios.



Real-time alarm with accurate location map

- Zone level
- Room level
- Floor level
- Indoor & outdoor

When an alarm is raised, DURAlarm displays the device location and other alarm details with both textual and graphical components. We can deploy the location accuracy, depending on client requirements and satisfying a mixture of duress needs across the campus.



Room Level Accuracy: Based on RTLS technology
The pop-up map when a duress event occurs in Room 6



Ward Level Accuracy: Based on the basestation location
The pop-up map when a duress event occurs in Ward 1

Open a voice channel allowing monitoring staff to listen in to the incident

DURAlarm can work with ENSUREspeech to open the voice channel enabling the response team to coordinate the appropriate response. All records will be stored for the use of reporting and training purposes.



Ascom Australia
2C Ground Floor, Building 2,
41-43 Bourke Rd
Alexandria NSW 2015.
Phone: + 02 9698 9000
ascom.com.au

