

# DATA SHEET



## Port Server (NUPS2/4-H)

### Features

- Combined module of a Room Gateway with medical input device
- Two SafeConnect sockets connect to medical devices via medical device smart cables
- Extendable up to four SafeConnect sockets
- Monitors intentional disconnects from the SafeConnect plug
- SafeConnect magnets use reverse polarity compatible with medical device smart cables
- SafeConnect blue background indicates medical cable connectivity only
- Easy mounting on industry standard VESA mount 75 x 75 mm
- Easy mounting on medical rail DIN 10 x 25 mm
- Easy mounting on medical rail Scandinavian Standard 10 x 30 mm
- Single RJ45 Socket for PoE+ interface
- Powered over Ethernet (PoE+ 802.3af)
- Easy to clean

## Technical Specifications

### NUPS2/4-H Overview

Item	NUPS2/4-H
Material	PBT/PC, ASA/PC enclosure rated UL94-V0
Color	Housing: White (NCS S0300-N / RAL9002)
Electrical connections	1x RJ45 socket for PoE+ interface 2x Magnetic Ascom SafeConnect receptacles (35N breakaway force) Optionally extendable to 4x Magnetic Ascom SafeConnect receptacles
LAN interface	100 BaseT Ethernet (auto negotiating)
Power	Input: PoE+ 802.3af <b>Note:</b> For installations requiring EMC Class B compliance, a PoE+ switch rated Class B must be used. For medical installations within 1.5 m of the patient, a PoE+ switch compatible with medical installations must be used (refer to installation guide for further information). Usage: max. 11.1 W
Wire type	Cat 5/5e/6/7, U/UTP
Isolation	1x MOPP (Means of Patient Protection)
Optical feedback	Red/Green LED
Buzzer	2400Hz, 65dB at 2.0m (6.6ft.)
Dimensions (W x H x D)	210 x 124 x 57 mm (8.27 x 4.88 x 2.24 in) (without medical rail clamp)
Weight	500 g (17.6 oz) (without medical rail clamp) 500 g (17.6 oz) medical rail clamp only
Mounting	Industry VESA 75 x 75 mm Medical Rail DIN 10 x 25 mm Medical Rail 10 x 30 mm

### Compatibility

Software Compatibility	Telligence 7.0 and later
------------------------	--------------------------

### Environmental

Operating temperature:	0°C to 50°C (32°F to 122°F)
Storage temperature:	-25°C to 55°C (-13°F to 131°F)
Relative Humidity:	30 to 95% (non condensing)

*Specifications are subject to change without notice*

Enclosure protection rating	IP40
-----------------------------	------

### Accessories

Accessories:	Medical Device Input Module (P/N NUMDI2-HE)
	SafeConnect Medical Device Smart Cable (P/N NUMDSC-H)
	Blank PD (P/N NUBPD-HE)
	Dummy Plug (P/N NU DP-H)
	Medical rail clamp DIN 10 x 25 mm (P/N NUMRC10x25-H)
	Medical rail clamp Scandinavian Standard 10 x 30 mm (P/N NUMRC10x30-H)

*Specifications are subject to change without notice*

## Compliance to Regulations and Standards

### Declaration of Conformity

The following complies with the relevant dates and amendments of all the following standards (valid references to dates and amendments are stated in the Declaration of Conformity):

### Regulatory Compliance EU/EFTA



EN 62368-1  
 EN 63000 (RoHS)  
 REACH  
 EN 60601-1  
 EN 60601-1-2 Class A

### Regulatory Compliance UK



Basic compliance of EU/EFTA standards if not superseded by UK regulations.  
 UK RoHS  
 HTM 08-03  
 EN 60601-1  
 EN 60601-1-2 Class A

### Regulatory Compliance ANZ



AS/NZS 62368-1  
 EN 60601-1  
 EN 60601-1-2 Class A

*Specifications are subject to change without notice*

**Ascom (Sweden) AB**

Grimbodalen 2

SE-417 49 Göteborg

Sweden

Phone +46 31 55 93 00

[www.ascom.com](http://www.ascom.com)

**ascom**